Instruction for Use of Sugita Titanium Clip II Applier

EN

Important information

The device is supplied non-sterile.


It must be reprocessed (cleaned and disinfected), using the validated reprocessing conditions. In addition, prior to use the device must be sterilized according to the validated sterilization conditions.

Specifications

Sugita Titanium Clip II Applier is exclusively designed for Sugita Titanium Aneurysm Clip II.

Sugita Titanium Clip II Appliers are intended as instruments used to hold and apply the Sugita Titanium Aneurysm Clip II during treatment of cerebral aneurysms.

"Sugita Titanium Clip II Applier" is hereinafter referred to as "applier".

[Figure 1]

Table 1. Product list

Code No.	Product Description
17-012-51G	Sugita Titanium Clip II Applier, Straight 70mm for Standard Clips
17-012-52G	Sugita Titanium Clip II Applier, Straight 90mm for Standard Clips
17-012-53G	Sugita Titanium Clip II Applier, Straight 110mm for Standard Clips
17-012-54G	Sugita Titanium Clip II Applier, 15°Angled for Standard Clips
17-012-55G	Sugita Titanium Clip II Applier, 30°Angled for Standard Clips
17-012-58G	Sugita Titanium Clip II Applier, Sideward Angled for Standard Clips
17-013-51G	Sugita Titanium Clip II Applier, Straight 70mm for Mini Clips
17-013-52G	Sugita Titanium Clip II Applier, Straight 90mm for Mini Clips
17-013-53G	Sugita Titanium Clip II Applier, Straight 110mm for Mini Clips
17-013-54G	Sugita Titanium Clip II Applier, 15°Angled for Mini Clips
17-013-55G	Sugita Titanium Clip II Applier, 30°Angled for Mini Clips

17-013-58G	Sugita Titanium Clip II Applier, Sideward Angled for Mini Clips
17-014-52G	Sugita Titanium Clip II Applier, Straight 90mm for Long Clips
17-014-54G	Sugita Titanium Clip II Applier, 15°Angled for Long Clips

Material: Titanium alloy, stainless steel

Types

Three types of appliers are available for Standard, Mini and Long clip each.

Appliers are color-coded on the applicable clips as follows. [Figure 1]

○ Appliers for Standard type: Distal part of appliers colored green.

○ Appliers for Mini type: Distal part of appliers: colored pink.

○ Appliers for Long type: Distal part of appliers: colored blue.

Applicable clip types are imprinted on the return spring of the appliers for your verification before use. [Figure 2]

e.g. "FOR SUGITA Ti-CLIP-II Standard ONLY"

[Figure 2]

Mechanism of action

Sugita Titanium Clip II Appliers operate under the following principles:

- When the Handle is opened and closed, Jaw opens and closes.
- Squeeze the Applier handles and set a clip between the Applier jaws by placing Sugita clips in the holding recesses of the jaws. When the handles are further squeezed with the clip in the jaws, the ratchet will be released to open the clip blade to the maximum width. At this point, gradual loosening of the handles will open the jaws and release the clip while closing the clip blades.

Intended purpose

"Sugita Titanium Clip II Applier" is an instrument used to hold and apply the Sugita Titanium Aneurysm Clip II for treatment of cerebral aneurysm.

Indications

The Sugita Titanium Clip II Applier is indicated for use in surgery where cerebral aneurysm clips are used.

Intended Users

The applier is intended for use by neurosurgeons who are thoroughly familiar with the surgical technique and the use of aneurysm clips and their appliers.

Intended patient population

The target treated population are adult patients, ≥18, for whom craniotomy clipping of cerebral aneurysms is indicated.

The performance characteristics of the device

Successfully apply the Sugita Titanium Aneurysm Clip II

Clinical benefit

Facilitate the safe placement of the Sugita Titanium Aneurysm Clip II

Directions for Use

Before using this product, inspect, wash, and sterilize in accordance with these instructions.

- 1) Combination clip with appliers
"Sugita Titanium Aneurysm Clip II"

Table 2. Sugita Titanium Aneurysm Clip II Color

Permanent Type		Temporary Type	
Type	Color	Type	Color
Standard Type	Spring - Green	Standard Type	Spring - Green
	Blade - Silver		Blade - Yellow
Mini Type	Spring - Pink	Mini Type	Spring - Pink
	Blade - Silver		Blade - Yellow
Long Type	Spring - Blue		
	Blade - Silver		

The clip is only used by the applicable clip applicator. Pay attention to the type of clip.

The applicable clip has been displayed on the clip applicator. Therefore, please use the appropriate clip.

- 2) By squeezing the applicator handles, the ratchet may be locked. Further squeezing of the handles will release the locked ratchet.
- 3) Set a clip between the applicator jaws by placing the clip correctly in the holding recesses of the jaws while the ratchet is released, then lock the ratchet. The clip will be held firmly between the jaws.
- 4) When the handles are further squeezed with the clip in the jaws, the ratchet will be released to open the clip blade to the maximum width. At this point, gradual loosening of the handles will close the clip blades with opening of the jaws to release the clip.
- 5) When holding the clip, it must be verified that the spring of the clip is correctly placed in the holding jaws of the clip applicator as shown in [Figure 3-A]. Improper holding of the clip may cause damage or ejection from the clip applicator. As this may present a risk during surgery, ensure correct holding of the clip. [Figure 3-B/C/D/E]
- 6) Do not hold the clip in angled positions.
- 7) Sugita Titanium Clip II Applicator is exclusively designed and constructed for use with Sugita Titanium Aneurysm Clips II only. It must NOT be used for clips of other manufacturers or conventional cobalt alloy clips.

A. Correctly positioned

B. Incorrectly positioned.

C. Incorrectly positioned

D. Incorrectly positioned

E. Incorrectly positioned

[Figure 3]

Contraindications

- 1) This applicator must not be used in patients with serious allergies to metals or exogenous materials.
- 2) This applicator must not be used in patients with infectious diseases or sepsis as a complication.

Warning

- 1) If the device is used in a patient with or is suspected of having CJD the instrument should be disposed of in accordance with local requirements.
- 2) Handling of this clip:
Rough handling of this applicator may cause deformation or damage to the applicator and may harm the patient.

- 3) Prohibition of use for unintended purposes

The applicator should not be used for other than intended purpose.

- 4) Clip applicators should not be used in MR environment.
- 5) If needed, all metal products or fragments thereof can be located by means of an X-Ray.

Caution

- 1) Use of this applicator is restricted to intended users only.
- 2) Do not use clips manufactured by other companies or incompatible with the applicator.
[Using the wrong clips may damage the aneurysm clip.]
- 3) Modifying or adulteration of this applicator is not allowed.
- 4) The applicator must be sterilized by steam sterilization, other methods should not be used.
- 5) Do not modify or reshape the return spring.
Appropriate width of the opening at holding clip is ruined by breaking and transforming the return spring and the clip might come to interfere, and eventually not come off easily with applicators.

Defects and Adverse Events (residual risk)

- 1) Serious adverse events
 - Infection
 - Damage to blood vessel
 - Death
 - Disease transmission
 - Cranial nerve damage
 - Blood loss
 - Organ damage
 - Aneurysm failure or aneurysm rupture.
 - Biological toxicity symptoms.
 - Biological inflammatory response
- 2) Other adverse events
 - General surgical complications
 - Delay in surgery

Notice

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in which the user and/or patient is established.

Summary of safety and clinical performance

For summary of the clinical performance and safety, see Summary of Safety and Clinical Performance (SSCP).

<https://ec.europa.eu/tools/eudamed>

Storage – Care and Handling

Store any new or unused products in a dry, clean, and safe place.

Validated Cleaning and Sterilization Procedure

The applicator is supplied in a non-sterile state.

It MUST be sterilized by the user (medical professional) prior to use using validated sterilization processes in accordance with guidelines stipulated by the government.

MIZUHO have confirmed safety and performance of the applicator up to 500 times of repeated cleaning and sterilization processing.

Proper sterilization and handling of this applicator is important for assuring that the applicators meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Sterilizers must be validated to ensure proper sterilization. Users (medical professionals) are responsible for validating the sterilization process and ensuring sterility of the individual sterilizer.

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer. Other procedures for minimizing contamination in the operating room, including the use of a surgical unit with laminar flow equipment at the

doctor's discretion, are also important.

Be sure to clean the applier prior to sterilization.

Bioburden prior to sterilization must be reduced to ensure successful sterilization.

Do not use chlorine-based solutions when cleaning this applier. This could cause metallic corrosion.

1) Cleaning

Important notes

- Clean the device immediately after use.
- Do not contact with other instruments.
- Use a detergent appropriate for this product and in accordance with the manufacturer's instructions.
- Cleaning is in two stages. Precleaning followed by automatic cleaning.

Conditions of cleaning validated by the manufacturer are as follows.

a) Precleaning

- ① Soak the product in detergent*^b at room temperature for 15 minutes, covering all surfaces. The temperature during cleaning should not exceed 40 °C. (20 °C to 40 °C)
- ② After the soaking time, this product is brushed for 1 minute in the solution. (while opening and closing the JAW part.)
- ③ Finally, this product is rinsed for 1 minute with running tap water.

b) Automated cleaning, disinfection and drying (washer disinfectant)

- Use the washer-disinfectant that has been officially recognized for its performance (according to the ISO 15883 series, with the CE mark, or approved by the FDA), install it accurately, check the operating status. Perform regular maintenance and inspection.
- After manual cleaning, perform cleaning according to the conditions in the table below.
- The ratchet should not be locked. (Conducted with release the ratchet)
- Make sure that the devices do not come into contact with each other.

Table 3. Washer-disinfectant recommended (validated) conditions

Stage	Procedure	Temp (°C)	Time (min)	Quality of water	Remarks
I	Preclean	<25	3	Tap water (Drinking water)	-
II	Washing	55	10	Deionised water * ^a	Detergent * ^b
III	Rinsing	<25	1	Deionised water * ^a	-
IV	Rinsing	<25	1	Deionised water * ^a	-
V	Hot water disinfection	90	5	Deionised water * ^a	-
VI	Drying	-	-	-	Follow the device program

*^a: Deionized water

Sufficient consideration must be given to the quality of the water used for diluting cleaning agents and the rinsing water used for medical device. When using pure water, highly purified water or sterile water for rinsing, we recommend rinsing water with a maximum viable count of 10 CFU/ml and endotoxin concentration of 0.25 EU/ml or less. If hard water contains metallic soap, or if the concentration of contaminant microorganisms or endotoxin is high, it may cause staining of the device and impair the cleaning and disinfection effect. (CFU: Colony Forming Unit)

*^b: Cleaning agent information used in verification.

Precleaning: Dr. Weigert, neodisher® MediClean forte, 0.5 %

Automated cleaning: Dr. Weigert, neodisher® MediClean forte, 0.4 %

- Use only slightly alkaline, low-foam detergent (detergent for aluminium, etc.) that will not impact the anodizing treatment of titanium.

2) Inspection

2)-1 Preinspection actions

If the product is wet or damp, dry it.

Use a water-based anticorrosive lubricant for medical devices.

MIZUHO recommend neodisher® IP Spray Dr. Weigert, REF: 4304 90.

Lubricating oil is completely removed by cleaning. After cleaning this product, apply a water-based, anticorrosive lubricant prior to sterilization. Do not use this product without applying anticorrosive lubricant on its sliding part. [Galling could occur.]

2)-2 Do not use the device if any abnormalities are observed by visual inspection or motion check.

- Damage including any sharp surfaces
- Deformation
- Corrosion
- Remaining dirt
- Corrosion
- Loosening of joints and screws
- Backlash of movable part
- Target Sugita clips cannot be properly hold.

Since the color of the surface may be faded in repetitive processing, do not use appliers that cannot be clearly distinguished with color coding. Please contact your distributor for repairs.

3) Packaging

This product can be sterilized by sterilizing the applier alone or by placing it in a dedicated sterilization case.

When using a sterilization case, store the product in the dedicated MIZUHO sterilization case (Table 4). In either case, make sure the ratchet has been released. Otherwise the return spring will be damaged.

When sterilizing the applier alone or in the dedicated sterilization case, double wrap the device with wraps conforming to ISO 11607 that are appropriate for steam sterilization.

Table 4. Sterilization case list

Code No.	Product Description
17-010-80	Sugita Titanium Clip II Case A for Appliers
17-010-06	Sugita Titanium Clips II Half Case
17-010-20	Sugita Titanium Clip II Metal Case A for Appliers
17-012-26	Sugita Titanium Clip II Metal Half Case

4) Sterilization

High-pressure steam sterilization is recommended for sterilization. The conditions in the table below have been validated.

Other sterilization methods are not recommended because the surface fades and affects the properties.

The high-pressure steam sterilizer should be verified, maintained and inspected according to ISO17665-1.

The ratchet should not be locked. (Conducted with release the ratchet)

Table 5. Validated Sterilization Method

Standard followed	ISO 17665-1	
Sterilization method	Pre-vacuum steam sterilization (Autoclave sterilization)	
Sterilization conditions	Sterilization temp.	Retention time
	132°C	4 minutes or more
	134°C	3 minutes or more
Drying conditions * ^a	Drying time (set time): 20 minutes (8 times fractionated vacuum)	

*^a: The actual drying time required depends directly on the parameters that are the sole responsibility of the user (e.g. load configuration and density, sterilizer conditions).

The high-pressure steam sterilizer used for verification is EN285-compliant.

5) Storage

Following sterilization, store in a clean place free from dust and moisture. Avoid high temperature and humidity.

Warranty

MIZUHO Corporation will repair defective parts of this product without charge for one year from the date of delivery/installment except for cases of damage caused by a third party's repair, act of nature, improper use or intentional damage. All other warranty terms and conditions are subject to regulations of MIZUHO Corporation.

Disposal

If the device needs to be disposed of, e.g. having been used on a patient with or suspected of having CJD, disposal shall be conducted in accordance with local regulations.

MIZUHO Corporation
3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan
<https://www.mizuho.co.jp>

EMERGO EUROPE
Westervoortsedijk 60, 6827 AT Arnhem
The Netherlands

Explanation/translation of symbols

	Catalogue number
	Batch code
	Manufacturer
	Date of manufacture
	Authorized representative in the European Community/European Union
	Importer
	Medical device
	Consult instructions for use or consult electronic instructions for use
	Keep dry
	Keep away from sunlight
	Non-sterile
	MR Unsafe
	Unique device identifier
	Do not use if package is damaged and consult instructions for use

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