

Be sure to read this instruction manual before using this product.  
Keep this manual available for reference when needed.

## Instructions for Use

Trade Name: Meter for Tourniquet

**Warning**

- Please read these instructions carefully prior to using this product. Product should be used according to these instructions, and pay close attention to the safety of patients. Not doing so may give rise to serious problems or adverse events.

**For the US market**

Do not reuse the device when it is used on a patient with or suspected of having Creutzfeldt-Jakob Disease (CJD) or variant CJD (vCJD).


**For the market outside the US**

When the device is used on a patient with or suspected of having CJD or vCJD, make sure to adhere to the most recent and updated restrictions available in each country and/or region for its reuse. Be cautious of the possibility of secondary infection. Refer to [www.a-k-i.org](http://www.a-k-i.org) or AAMI Standard ST79 for more information related to cleaning and sterilization.

**Contraindication / Prohibition**

- Use for intended purpose only  
Use devices for their intended purposes only. Incorrect use could cause this product to break.
- Use with specified products only  
Use this product only with products specified by Manufacturer. Other products than those specified by Manufacturer could be incompatible with this product due to differences in design and development policies.
- Prohibition of use of chemicals  
Avoid exposing this product to chemicals. Doing so could damage this product due to corrosion.
- Handle with care  
Handle this product with care, as it can be deformed or damaged. Rough handling could significantly reduce the service life of devices and appliances.
- Do not sterilize  
Do not sterilize this product. [Doing so may affect the capabilities of the product.]

**Shape / Structure**



Item	Code No.	Product Description
A	01-440-06	Meter for Tourniquet
B	01-447-25	Tourniquet Joint Hose C

	Composition	Material
Meter for Tourniquet Code No.:01-440-06	Main body	Stainless Steel
		Aluminum alloy
		Glass
		Brass
		Iron
		Rubber
Tourniquet Joint Hose C Code No.:01-447-25		Stainless Steel
		Resin

**Intended Purpose**

A band like device for reusable purposes which encircles an extremity (arm or leg) to control the circulation and prevent normal blood flow to or from the distal area. It is used in combination with a tourniquet unit which controls the pressure applied. The cuff usually comprises a dual-chamber allowing for alternation of the pressure site.

**Instructions for use**

-Before using this product, inspect, wash, and sterilize in accordance with these instructions

**CAUTIONS**

- Defect/Adverse event  
Defect
  - Corrosion and pitting caused by use of chemicals
  - Damage or breakage caused by the corrosion or pitting
- Pressure decrease  
Pressure may decrease due to a change in the patient's posture. In such cases, apply pressure again.
- Allergic symptoms  
There is no direct contact, but the meter includes natural rubber parts. In rare circumstances, natural rubber may cause allergic reactions such as itching, redness, hives, and swelling. If symptoms occur, terminate use immediately and perform proper treatment.

**Storage / Life**

- Please store the device in normal ambient temperature areas. Do not store in areas of high humidity where the temperature may dramatically vary causing condensation. Do not store on or near chemicals as the chemicals may cause damage to the device.
- Service life of this product: 5 years  
(Subject to following manufacturer's specified maintenance, inspection and proper storage requirements.)

**Maintenance / Inspection**

- Check prior to each use  
Operational and functional checks  
Conduct daily and pre-operation checks of this product to make sure that it functions properly.
- Check after each use
  - Cleaning and disinfecting  
Remove blood, medicine, or dirt with water, and disinfect product with alcohol or gauze.

**Packing**

1 piece per pack

**Maintenance and check by agents**

For safe use of this instrument, inspect device prior to and after each use. Alternative or no review and/or inspection may cause injury to the patient and/or healthcare worker and may decrease the performance and function of this device. Additionally, it is recommended to schedule a periodic inspection through your local distributor or the manufacturer.

**Matters related to warranty period**

MIZUHO Corporation will repair defective parts of this product without charge for one year from the date of delivery/installment except for cases of damage caused by a third party's repair, act of nature, improper use or intentional damage. All other warranty terms and conditions are subject to regulations of MIZUHO Corporation.

**Name and address of manufacturer**

MIZUHO Corporation  
3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, Japan  
<http://www.mizuho.co.jp>