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INSTRUCTIONS FOR THE USE OF SUGITA AVM MICROCLIPS

Index

- | | |
|---------------------------------------|-------------------------------------|
| 1. Specifications | 10. Warning |
| 2. Types | 11. Caution |
| 3. Spring Type | 12. Defects and Adverse Events |
| 4. Action Mechanism | 13. Storage - Care and Handling |
| 5. Surface Finish | 14. Sterilization |
| 6. Blade Cross-sectional Shape | - Wrapping Prior to Sterilization - |
| 7. Indication for Use | - Sterilization - |
| 8. Directions for Use | - Storage - |
| 9. Contraindications and Prohibitions | |

1. Specifications

Sugita AVM Microclips are made of cobalt-chromium-nickel-molybdenum-iron alloy (ISO 5832-7/ASTM F1058 G2). "Sugita AVM Microclip" herein referred to as "clip".


ASTM: American Society for Testing and Materials

ISO: International Organization for Standardization

This clip designed for **single-use only** is exclusively indicated only for use with the Sugita AVM Microclip Applying Forceps (hereinafter this is called "clip applier"), and requires sterilization prior to use. Reuse may cause breakage of clips and may cause injury or infection.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a measuring point located 1.0 mm from the distal end of the blade length (L) of the clip [Figure 1] while its blades are opened at a 1.0 mm width. The result of this measurement is shown in newtons and grams on the package label.

The product information including name and size are as indicated on the clip's package label.

[Figure 1]

2. Types

Prod. Code	Length of blade
07-932-02	2 mm
07-932-03	3 mm
07-932-04	4 mm
07-932-05	5 mm

For specific dimensions and shapes, please see this clip's catalog.

3. Spring Type : Integral Single Bend

4. Action Mechanism : Alligator Action

5. Surface Finish

- Clip Body: Satin
- Material Treatment: Electrolytic Polish

6. Blade Cross-sectional Shape : Rolled

7. Indication for Use

This device is indicated for temporary occlusion of cerebral blood vessels.

For the US market

This device is indicated for temporary occlusion of blood vessels.


8. Directions for Use

This clip is packaged clean, non-sterile, and must be properly sterilized before use.

Each clip must be held with the specified clip applier.

When holding the clip, it must be verified that the spring of the clip is correctly placed in the holding jaws of the clip applier as shown in [Figure 2]. Improper holding of the clip may cause damage or ejection from the clip

applier. As this may present a risk during surgery, ensure correct holding of the clip.

[Figure 2]

- Do not hold the clip in angled positions.
- If the spring of the clip is not properly held with the holding jaws of the clip applier, the clip may slip or spring out from the jaws. This may also damage the spring of the clip, and there is a risk of breakage and reduced holding force, which may prevent closure of the targeted cerebral blood vessels.

It is recommended that doctors record the type, lot, and serial number of the clip used in the operation record for each patient.

To support traceability, Patient/Trace labels are included in the clip package.

9. Contraindications and Prohibitions

1) Prohibition of reuse

This clip is for single use only. The clip must be discarded once the blades of the clip have been opened or the clip has come into contact with the blood or tissue of a patient. Reuse may cause breakage of clips and may cause injury or infection.

2) Prohibition of use for unintended purposes

This clip should not be used for purposes other than indication for use.

The clip should not be used for permanent occlusion of, or permanent placement in blood vessels.

3) To prevent changes in the clips closing force or in the mechanical features of the clip when it is handled, only use the exclusive clip applier. Do not use clip applicators from other manufacturers.

4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when spring of the clip is placed in the proper position inside the clip applier.

5) Modifying or adulteration of this clip is not allowed.

6) In order to avoid a galvanic current, do not use this clip with implants from other manufacturers or with implants made of different materials.

A galvanic current may cause corrosion and weaken the holding force.

7) This clip must not be used in patients with serious allergies to metals or exogenous materials.

8) This clip must not be used in patients with infectious diseases or sepsis

-
- as a complication.
- 9) This clip must not be sterilized using low-temperature hydrogen peroxide gas plasma.

This clip is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma. [Sterilization using such method may cause surface color deterioration and affect its characteristics.]

10. Warning

- 1) It is recommended the patient is informed about the procedure and characteristics of this clip.
- 2) Handling of this clip:

Rough handling of this clip may cause deformation or damage to the clip and may harm the patient.

11. Caution

Use of this clip is restricted to surgeons with proper training in the use of this clip.

Impact should not be applied to its package or aluminum container in order to avoid clips falling out of its packaging.

- The aluminum container is sealed to maintain the clip integrity. Clips removed from their original packaging or sterilized by user cannot be returned.
 - This clip can be re-sterilized only if the sterilization bag or wrap is not opened in the operation room.
-

-
- If the sterilization bag or wrap is opened in the operating room, the clip must be discarded.
 - Provide a sterile barrier (pouch or wrap) to protect the clip and its container from contamination. Once the sterile barrier has been removed and the clip and its container has been subjected to possible contamination, the clip and its container must be disposed.

MR Conditional

Non-clinical testing demonstrated that this clip is MR Conditional. A patient with this clip can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only
- Maximum spatial gradient magnetic field of 1,000 gauss/cm (10-T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, this clip is expected to produce a maximum temperature rise of 2.4°C (4.32°F) after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by this clip extends

approximately 25-mm from this clip when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

NOTE: Do not take the clip applier into the MR environment.

The clip applier contains material which has been determined to be MR Unsafe.

12. Defects and Adverse Events

1) Serious adverse events

- Displacement of the clip (slippage)
- Ejection of the clip from clip applier
- Breakage of the clip
- Cerebrovascular injury
- Hemorrhage
- Death

2) Other adverse events

- Infection of surgical wound
- General surgical complications
- Delayed clip displacement due to the application of magnetic field
- Occlusion of adjacent blood vessel
- Foreign body reaction

13. Storage – Care and Handling

This clip is prepared in a non-sterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applier other than the exclusive clip applier. Handling of this clip should be minimized to reduce contamination or damage.

14. Sterilization

This clip is supplied non-sterile. Therefore, this clip must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

Proper sterilization and handling of this clip is important for assuring that the clips meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures for minimizing contamination in the operating room, including the use of a surgical unit with laminar flow equipment at the doctor's discretion, are also important.

Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.

Clean and sterilize the clip applier according to the clip applier instruction manual.

- Wrapping Prior to Sterilization -

Do not wash the clip. This clip is provided clean and the clip and its packaging is not designed to be washed or cleaned prior to sterilization by user.

Remove the sealed aluminum container from the outer box/original packaging and wrap the container with sterilization packaging material (nonwoven fabric or sterilization bag/pouch).

Place in an appropriate sterilization bag/pouch or wrap and sterilize. Do not remove the clip from the aluminum container.

To support sterilization, each individual clip is retained in a heat-resistant autoclavable foamed silicone material within an aluminum container.

- Sterilization -

〈Validated Sterilization Method〉

Standard followed	ISO 17665-1		
Sterilization method	Pre-vacuum steam sterilization (Autoclave sterilization)		
Sterilization conditions	Sterilization temp.	Retention time	Drying time
	132°C / 270°F*	4 minutes or more	20 minutes or more
	134°C / 273°F*	3 minutes or more	20 minutes or more

*According to the Guideline for Disinfection and Sterilization in Healthcare

Facilities, 2008 (Update: May 2019)

Do not sterilize using low-temperature hydrogen peroxide gas plasma. It may cause discoloration and/or deformation.

With a high-pressure steam sterilizer that has been confirmed to function properly, the validation of sterilization is verified by a medical institution, sterilization conditions with a proven effectiveness shall be applied.

- Storage -

After sterilization, the clip should be stored and protected in a clean and dry area.

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