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# INSTRUCTIONS FOR THE USE OF SUGITA AVM MICROCLIPS II

## Index

- |                                       |  |
|---------------------------------------|--|
| 1. Specifications                     | 10. Warning  |
| 2. Types                              | 11. Use of MRI                                     |
| 3. Spring Type                        | 12. Caution  |
| 4. Action Mechanism                   | 13. Defects and Adverse Events                     |
| 5. Surface Finish                     | 14. Validated Cleaning and Sterilization Procedure |
| 6. Blade Cross-sectional Shape        | - Sterilization -                                  |
| 7. Indication for Use                 | - Storage -  |
| 8. Directions for Use                 |  |
| 9. Contraindications and Prohibitions |  |

---

## 1. Specifications


Sugita AVM Microclips II are made of cobalt-chromium-nickel-molybdenum-iron alloy (ISO 5832-7/ASTM F1058 G2). “Sugita AVM Microclip II” herein referred to as “clip”.

This clip designed for **single-use only** is exclusively indicated only for use with the Sugita AVM Microclip II Applying Forceps (hereinafter this is called “clip applier”), and requires sterilization prior to use. Reuse may cause breakage of clips and may cause injury or infection.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a measuring point located 1.0 mm from the distal end of the blade length (L) of the clip [Figure 1] while its blades are opened at a 1.0 mm width. The result of this measurement is shown in newtons and grams on the package label.

The product information including name and size are as indicated on the clip's package label.

The clip is set in a dedicated individual case.

[Figure 1]

## 2. Types

Prod. Code	Length of blade
07-932-12	2 mm
07-932-13	3 mm
07-932-14	4 mm
07-932-15	5 mm

---

**3. Spring Type** : Integral Single Bend

**4. Action Mechanism** : Alligator Action

**5. Surface Finish**

- Clip Body: Satin
- Material Treatment: Electrolytic Polish

**6. Blade Cross-sectional Shape** : Rolled

**7. Indication for Use**

This clip is used to block blood flow in cerebral blood vessels.

**8. Directions for Use**


**Clips are provided in non-sterile packaging and should be properly cleaned and sterilized prior to use. Remove from the individual case using clip applicator just prior to clipping.**

**Each clip must be held with the specified clip applicator. Each clip applicator is etched to read; e.g. “for AVM Microclip II ONLY”.**

When holding the clip, it must be verified that the spring of the clip is correctly placed in the holding jaws of the clip applicator as shown in **[Figure 2]**. Improper holding of the clip may cause damage or ejection from the clip applicator. As this may present a risk during surgery, ensure correct holding of the clip.

---

〈Correct holding method〉

〈Improper holding method〉

a) The holding of the clip is shallow.

b) The holding of the clip is deep.

c) Angular misalignment

d) Angular misalignment

[Figure 2]

Following completion of the surgery, physicians must record the type, lot and serial number of each clip used in the operation record for each patient.

Patient record labels, Implant card are included in the clip package to

---


---

facilitate traceability.

Be sure to fill in the necessary information on the implant card and provide it to the patient.


〈How to remove the clip〉

• Individual case and clip set state **[Figure 3-1·2·3]**


[Figure 3-1]

- 
- Since the clip is fixed by sandwiching it between plates of the same material as the clip set in the inner case, the clip can be directly held with clip applicator and removed. [Figure 3-2•3]

[Figure 3-2]

- Check the clip set in the inner case and correctly pinch it in the groove of the holding jaws of the clip applicator to remove it. [Figure 3-3]

Hold the clip after applying the tip of the clip applicator to the inner case.

[Figure 3-3]

---

## 9. Contraindications and Prohibitions

### 1) Prohibition of reuse

This clip is for single use only. The clip must be discarded once the blades of the clip have been opened or the clip has come into contact with the blood or tissue of a patient. Reuse may cause breakage of clips and may cause injury or infection.

### 2) Prohibition of use for unintended purposes

This clip should not be used for purposes other than indication for use.

### 3) To prevent changes in the clips closing force or in the mechanical features of the clip when it is handled, only use the exclusive clip applicator. Do not use clip applicators from other manufacturers.

### 4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the proper holding force. The blades of clip should be opened only when spring of the clip is placed in the proper position inside the clip applicator.

### 5) Modifying or adulteration of this clip is not allowed.

### 6) In order to avoid a galvanic current, do not use this clip with implants from other manufacturers or with implants made of different materials. A galvanic current may cause corrosion and weaken the holding force.

### 7) This clip must not be used in patients with serious allergies to metals or exogenous materials.

### 8) This clip must not be used in patients with infectious diseases or sepsis as a complication.

### 9) This clip must not be sterilized using low-temperature hydrogen

---

---

peroxide gas plasma.

This clip is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma.

## 10. Warning

- 1) It is recommended the patient is informed about the procedure and characteristics of this clip.
- 2) Handling of this clip:  
Rough handling of this clip may cause deformation or damage to the clip and may harm the patient.

## 11. Use of MRI

MR Conditional

Non-clinical testing has demonstrated the “Sugita AVM Microclip II” are MR Conditional.

A patient with this device can be safely scanned in an MR system meeting the following conditions.

Failure to follow these conditions may result in injury.

---

<b>Parameter</b>	<b>Condition</b>
Device Name	Sugita AVM Microclip II
Static Magnetic Field Strength (B0)	1.5 T and 3 T
Maximum Spatial Field Gradient	16 T/m (1,600 gauss/cm)
RF Excitation	Circular Polarized (CP)
RF Transmit Coil Type	There are no Transmit Coil restrictions.
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Scan Duration	2 W/kg whole-body average SAR for 15 minutes of continuous RF (a sequence or back-to-back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact of 17.1 mm.

NOTE : Do not take the clip applicator into the MR environment.

The clip applicator contains material which has been determined to be MR Unsafe.

## 12. Caution

This product is designed for use by healthcare professionals, especially trained neurosurgeons with experience in its use.

## 13. Defects and Adverse Events

1) Serious defects

- Clip migration

- 
- Significant decrease in closing force
  - Misalignment of the blade tip
  - Inversion
  - Displacement of the clip (slippage)
  - Ejection of the clip from the clip applier
  - Breakage of the clip
- 2) Other defects
- Delayed clip displacement due to the application of magnetic field
- 3) Serious adverse events
- Rupture of a cerebral aneurysm
  - Cerebrovascular injury
  - Hemorrhage
  - Death
- 4) Other adverse events
- Infection of surgical wound
  - General surgical complications
  - Occlusion of adjacent blood vessel
  - Foreign body reaction
-

---

## **14. Validated Cleaning and Sterilization Procedure**

The clip is supplied in a non-sterile state. It should therefore be sterilized by the user (medical professional) prior to use using validated sterilization processes in accordance with guidelines stipulated by the government.

Clips are cleaned and sterilized while kept in individual cases. Do not open the individual cases during the cleaning and sterilization process.

Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.

Clean and sterilize the clip applicator according to the clip applicator instruction manual.

This clip and clip applicator must be cleaned to reduce bioburden before first sterilization and resterilization.

Do not use chlorine-based solutions when cleaning this clip. This could cause metallic corrosion.

### **1) Cleaning**

When cleaning, set the individual case to “07-932-11: SUGITA AVM Microclip II CASE B”.

---

⟨Washer disinfectant (automatic washing)⟩

Stage	Procedure	Temperature [°C (°F )]	Time [min]	Quality of water	Remarks
I	Preclean	<25 (77)	3	Tap water (drinking water)	–
II	Washing	55 (131)	10	deionised water <sup>*a</sup>	detergent <sup>*b</sup>
III	Rinsing	<25 (77)	1	deionised water <sup>*a</sup>	–
IV	Rinsing	<25 (77)	1	deionised water <sup>*a</sup>	–
V	Hot water disinfection	90 (194) (A <sub>0</sub> : 3000)	5	deionised water <sup>*a</sup>	–
VI	Drying	–	–	–	follow the device program

\*a: Deionized water

Sufficient consideration must be given to the quality of the water used for diluting cleaning agents and the rinsing water used for medical device. When using pure water, highly purified water or sterile water for rinsing, we recommend rinsing water with a maximum viable count of 10 CFU/ml and endotoxin concentration of 0.25 EU/ml or less. If hard water contains metallic soap, or if the concentration of contaminant microorganisms or endotoxin is high, it may cause staining of the device and impair the cleaning and disinfection effect. (CFU: Colony Forming Unit)

\*b: Cleaning agent information used in verification

Dr. Weigert, neodisher® MediClean forte, 0.4% (worst case)

---

## **2) Inspection**

Do not use if the following items are recognized in individual cases.

- Damage
- Deformation
- Remaining dirt

## **3) Packaging**

When sterilizing, wrap in one of the following.

- “Unopened individual case” alone.
- “Unopened individual case” + “07-932-11: SUGITA AVM Microclip II CASE B”
- “Unopened individual case” + “07-932-11: SUGITA AVM Microclip II CASE B” + “07-932-10: SUGITA AVM Microclip II CASE A”

For packaging, use the medical packaging material “non-woven fabric or sterilization bag” for high-pressure steam sterilization that has microbial barrier properties (conforms to ISO 11607).

---

## - Sterilization -

〈Validated Sterilization Method〉

Standard followed	ISO 17665-1		
Sterilization method	Pre-vacuum steam sterilization (Autoclave sterilization)		
Sterilization conditions	Sterilization temp.	Retention time	Drying time
	132°C / 270°F*	4 minutes or more	30 minutes or more
	134°C / 273°F*	3 minutes or more	30 minutes or more

\*According to the Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008 (Update: May 2019)

With a high-pressure steam sterilizer that has been confirmed to function properly, the validation of sterilization is verified by a medical institution, sterilization conditions with a proven effectiveness shall be applied.

## - Storage -

After sterilization, the clip should be stored and protected in a clean and dry area.

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
**4. 作動機構** : Alligator Action

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- Clip Body : Satin

- Material Treatment : Electrolytic Polish

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