SUGITA

Titanium ANEURYSM CLIP II BOOSTER

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INSTRUCTIONS FOR THE USE OF SUGITA TITANIUM ANEURYSM CLIP II BOOSTER

Index

- | | |
|---------------------------------------|---|
| 1. Specifications | 10. Warning |
| 2. Types | 11. Use of MRI |
| 3. Spring Type | 12. Caution |
| 4. Action Mechanism | 13. Defects and Adverse Events |
| 5. Surface Finish | 14. Storage - Care and Handling |
| 6. Blade Cross-sectional Shape | 15. Validated Cleaning and Sterilization
Procedure |
| 7. Indication for Use | - Validated Sterilization Procedure - |
| 8. Directions for Use | - Storage - |
| 9. Contraindications and Prohibitions | |

1. Specifications

The Sugita Titanium Aneurysm Clip II Booster is used to reinforce the closing force of Sugita Titanium Aneurysm Clips II Permanent Type (standard type) as shown in [Table 1].


“Sugita Titanium Aneurysm Clip II Booster” herein referred to as “clip”.

This clip has a channel on the inside of the blades. Attach this clip so that its channel is aligned with the blades of the clip whose closing force you want to increase and so that the springs are parallel [Figure 1 (A)]. Do not use if the springs are not in a parallel position [Fig. 1 (B)].

This clip should not be used with clips other than the standard Sugita Titanium Aneurysm Clips II Permanent Type.

[Table 1]

	Prod. Code
Applicable Type	17-001-** 01, 02, 05, 06, 09 to 11, 15, 17, 18, 20 to 25, 47 to 50, 60 to 63, 65, 66, 69, 70, 73 17-002-01 to 03

[Figure 1]

This device is made of titanium 6-aluminium 4-vanadium alloy (Ti-6Al-4V ELI; as specified by ISO5832-3/ASTM F136) and unalloyed titanium (specified by ISO 5832-2/ASTM F67 G2).


ASTM: American Society for Testing and Materials

ISO: International Organization for Standardization

This clip is designed for **single-use only**, is exclusively indicated only for use with the Sugita Titanium Aneurysm Clips II Applying Forceps (hereinafter this is called “clip applicer”), and requires cleaning and sterilization prior to use. Reuse may cause breakage of clips and may cause injury or infection.

In order to ensure the traceability, a lot number is shown on each clip.

The closing force of each clip is measured during the final quality assurance stage with a computerized electronic scale at a specified measuring point as shown in **[Figure 2]**. The result of this measurement is shown in newtons and grams on the package label.

[Figure 2]

2. Types

Permanent Type

Type	Color
Booster	Spring - Green
	Blade - Green

For specific dimensions and shapes, please see this

device's catalog.

3. Spring Type : Integral Coil

4. Action Mechanism : Alligator Action

5. Surface Finish

- Clip Body : Satin

- Material Treatment : Anodizing Treatment

6. Blade Cross-sectional Shape : Rolled

7. Indication for Use

This device is used to reinforce the closing force of standard Sugita Titanium Aneurysm Clips II for permanent occlusion (permanent clips).


8. Directions for Use

This device is packaged non-sterile, and must be properly sterilized before use.

This device must be held with the clip applicator.

When holding the clip, it must be verified that the spring of the clip is correctly placed in the holding jaws of the clip applicator as shown in [Figure 3].

Improper holding of the clip may cause damage or ejection from the clip applicator. As this may present a risk during surgery, ensure correct holding of the clip.

[Figure 3]

-
- Do not hold the clip in angled positions.
 - If the spring of the clip is not properly held with the holding jaws of the clip applicator, the clip may slip or spring out from the jaws. This may also damage the spring of the clip, and there is a risk of breakage and reduced holding force.

It is recommended that doctors record the type, lot, and serial number of the clip used in the operation record for each patient.

To support traceability, Patient/Trace labels are included in the clip package.

9. Contraindications and Prohibitions

1) Prohibition of reuse

This device is for single use only. The clip must be discarded once the blades of the clip have been

opened or the clip has come into contact with the blood or tissue of a patient. Reuse may cause breakage of the clip and injury or infection of the tissue.

2) Prohibition of use for unintended purposes

This device should not be used for purposes other than reinforcing the closing force of Sugita Titanium Aneurysm Clips II Permanent Type (standard type). Do not use this device with clips other than the standard Sugita Titanium Aneurysm Clips II Permanent Type.

3) To prevent changes in the clips closing force or in the mechanical features of the clip when it is handled, only use the exclusive clip applicator. Do not use clip applicators from other manufacturers.

4) Do not open the blades of clip directly with fingers. This may cause corrosion or alter the

proper holding force. The blades of clip should be opened only when spring of the clip is placed in the proper position inside the clip applicator.

- 5) Modifying or adulteration of this product is not allowed.
- 6) In order to avoid a galvanic current, do not use this device with implants from other manufacturers or with implants made of different materials.
A galvanic current may cause corrosion and weaken the holding force.
- 7) This device must not be used in patients with serious allergies to metals or exogenous materials.
- 8) This device must not be used in patients with infectious diseases or sepsis as a complication.
- 9) This device must not be sterilized using low-temperature hydrogen peroxide gas plasma.

This product is not compatible with sterilization using low-temperature hydrogen peroxide gas plasma. [Sterilization using such method may cause surface color deterioration and affect its characteristics.]

10. Warning

- 1) It is recommended the patient is informed about the procedure and characteristics of this clip.
- 2) Closing status of the blades of clip should be verified immediately following surgery.

[The most serious adverse events that may be reported in conjunction with this clip include displacement of the clip, slip-out of the clip, ejection of the clip from the clip applier, clip fracture, laceration of a cerebral aneurysm,

cerebrovascular injury, hemorrhage, death, etc.]

3) Handling of this product:

Rough handling of this product may cause deformation or damage to the clip and may harm the patient.

11. Use of MRI

MR Conditional

Non-clinical testing demonstrated that this device is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla, only

-
- Maximum spatial gradient magnetic field of 1,000 gauss/cm (10-T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, this device is expected to produce a maximum temperature rise of 2.2°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by this device extends approximately 5-mm from this device when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

NOTE: Do not take the clip applicator into the MR environment.

The clip applicator contains material which has been determined to be MR Unsafe.

12. Caution

Use of this product is restricted to surgeons with proper training in the use of this product.

13. Defects and Adverse Events

1) Serious adverse events

- Displacement of the clip (slippage)
- Ejection of the clip from clip applicator
- Breakage of the clip
- Rupture of a cerebral aneurysm
- Cerebrovascular injury
- Hemorrhage
- Death

2) Other adverse events

- Infection of surgical wound
- General surgical complications
- Delayed clip displacement due to the application of magnetic field
- Occlusion of adjacent blood vessel
- Foreign body reaction

14. Storage – Care and Handling

This device is prepared in a non-sterile package and should be stored in its package in a cool, dry, and secure location.

Do not use any clip applier other than the exclusive clip applier.

Handling of this device should be minimized to reduce contamination or damage.

15. Validated Cleaning and Sterilization Procedure

This device is supplied non-sterile. Therefore, this device must be sterilized by the users (healthcare professionals) prior to surgery using validated sterilization processes in accordance with guidelines stipulated by the government.

Proper sterilization and handling of this clip is important for assuring that the clips meet sterile specifications prior to surgery as an aid in minimizing the occurrence of infection in the surgical patient.

Other procedures for minimizing contamination in the operating room, including the use of a surgical unit with laminar flow equipment at the doctor's discretion, are also important.

Sterilization equipment should be validated to assure that it sterilizes specific load configurations properly. The end user is responsible to validate the sterilization process in the individual sterilizer and assure the sterility.

Sterilization equipment should be cleaned and maintained in accordance with standardized operating procedures specified by its manufacturer.

Clean and sterilize the clip applicator according to the clip applicator instruction manual.

This clip and instruments must be cleaned to reduce bioburden before first sterilization and re-sterilization.

Do not use chlorine-based solutions when cleaning this product. This could cause metallic corrosion.

1) Cleaning

Put the clip in the designated Sugita Clip Case or Sugita Clip Tray (hereinafter this is called “clip case or tray”) and close the lid.

<In case of using Ultrasonic cleaner>

(1) Fully immerse the case (tray) in the water (tap water) of the washing tank.

(2) Ultrasonic cleaning

Cleaning time: 15 minutes or more

Detergent: Neutral enzyme type detergent (Diluted concentration 1.0%)

(3) Rinsing

- Make sump water with running water, immerse the case (tray) completely in water (tap water or deionized water), shake the case (tray) 5 times up and down.

- Renew the water and shake the case (tray) 5 times up and down in the water again.
- Perform a series of rinsing operations three or more times in total.

(4) Dry thoroughly.

Recommended drying condition: 110 °C, 30 minutes or more

<In case of using washer disinfectant (automatic washing)>

Stage	Procedure	Temperature [°C]	Time [min]	Quality of water	Remarks
I	Preclean	<25	3	Tap water	–
II	Washing	55	10	Tap water	Alkaline detergent (such as cleaning agent for aluminum) that does not affect anodic oxidation coatings Diluted concentration 0.5%
III	Rinsing	<25	4	Deionized water	Set rinsing process program as “1 time (1 minute) × 4 times = 4 minutes”
IV	Hot water disinfection	–	–	Deionized water	Follow the program of the device.
V	Drying	–	–	–	

2) Inspection

Check products visually one by one and do not use if the following abnormalities are found.

- Damage
- Deformation
- Corruption
- Inappropriate blade position (deviation, clearance, tip unevenness, etc.)
- Remaining dirt

Since the color of the surface may be faded in repetitive processing, do not use products that cannot be clearly distinguished with color coding.

3) Packaging

Place clip in the designated clip case or tray and close the lid.

Then pack the case (tray) containing the product using standard medical packaging material (non-woven fabric or sterilized bag) with protective sterile barrier complying with ISO 11607 for steam sterilization.

- Validated Sterilization Procedure -

With the pre-vacuum type high-pressure steam sterilizer which functions properly, this product shall be sterilized so that sterilization conditions [134 °C for 3 minutes or more] are maintained.

Do not sterilize using low-temperature hydrogen peroxide gas plasma. It may cause discoloration and/or deformation.

With a high-pressure steam sterilizer that has been confirmed to function properly, the validation of

sterilization is verified by a medical institution, sterilization conditions with a proven effectiveness shall be applied.

- Storage -

Do not store the device in high temperatures, or in areas with high humidity and/or where the temperature varies. Please store the device in clean dust-free and dry area.

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