

Lithotomy Positioning System

Levigator Stirrup II



9 Reasons to Select Mizuho Stirrups

Interchangeable boot design saves storage space and cost

Quick attachment / release with integrated clamps

Easy and smooth adjustment with gas cylinder and ergonomic handle

Wide range of motion

Lithotomy: +80° to -35°
Abduction: +30° to -10°

Lightweight design: 6.5 kg per unit

*regular boot, pad and clamp included

Patient weight capacity: 270 kg

Contoured boots with soft and thick pads

Made-in-Japan quality

Compatible with JP, US, and EU side rails

Regular Stirrup Set

Main Unit + Boot + Pad

	Product Code	Description
R	1 8 - 0 7 0 - 9 1	Levigator Stirrup II Set, Regular, Right Leg
L	1 8 - 0 7 0 - 9 2	Levigator Stirrup II Set, Regular, Left Leg

*Right and left units are sold separately



Maximize positioning flexibility while minimizing storage space and investment

- A variety of boot sizes and shapes accommodate different patient types
- Interchangeable boots design saves storage space and cost

Main Unit



	Product Code	Description
R	18-070-93	Levigator Stirrup II Main Unit, Right Leg
L	18-070-94	Levigator Stirrup II Main Unit, Left Leg

*Sold separately

Levigator Cart




Product Code	Description
18-070-61	Levigator Stirrups Cart

Short Boot Extender



Provides additional length when needed, to help achieve proper foot positioning with the Short Boot



Product Code	Description
18-070-50	Short Boot Extender

*Order two pieces for both legs

Regular

Standard adult size



	Product Code	Description
R	18-070-21	Boot for Levigator Stirrup II, Regular, Right Leg
L	18-070-22	Boot for Levigator Stirrup II, Regular, Left Leg
R	18-070-41	Boot Pad for Levigator Stirrup II, Regular, Right Leg
L	18-070-42	Boot Pad for Levigator Stirrup II, Regular, Left Leg

*Sold separately

Fin Type

Lateral knee support is designed to help reduce the risk of nerve injury



	Product Code	Description
R	18-070-23	Boot for Levigator Stirrup II, Fin Type, Right Leg
L	18-070-24	Boot for Levigator Stirrup II, Fin Type, Left Leg
R	18-070-43	Boot Pad for Levigator Stirrup II, Fin Type, Right Leg
L	18-070-44	Boot Pad for Levigator Stirrup II, Fin Type, Left Leg

*Sold separately

Wide Type

For larger and bariatric patients



	Product Code	Description
R	18-070-27	Boot for Levigator Stirrup II, Wide, Right Leg
L	18-070-28	Boot for Levigator Stirrup II, Wide, Left Leg
R	18-070-47	Boot Pad for Levigator Stirrup II, Wide, Right Leg
L	18-070-48	Boot Pad for Levigator Stirrup II, Wide, Left Leg

*Sold separately

Short Type

For pediatric and small patients



	Product Code	Description
R	18-070-25	Boot for Levigator Stirrup II, Short, Right Leg
L	18-070-26	Boot for Levigator Stirrup II, Short, Left Leg
R	18-070-45	Boot Pad for Levigator Stirrup II, Short, Right Leg
L	18-070-46	Boot Pad for Levigator Stirrup II, Short, Left Leg

*Sold separately



MIZUHO Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN



Website



Contact Us



EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

* The information and specifications are subject to change without prior notice due to ongoing technical improvement.

* All CG images and illustrations are for reference only.

©MIZUHO 2026

2606-PDF-MC-EN/BR-001-2606