




OPERATING TABLE ACCESSORIES

C O N T E N T S

| | | |
|-----------|--|-----------|
| 1 | ARM SUPPORT | 1 |
| 2 | LEG SUPPORT | 9 |
| 3 | HEAD SUPPORT | 15 |
| 4 | BODY SUPPORT | 18 |
| 5 | X-RAY ACCESSORIES | 20 |
| 6 | ORTHOPEDICS / TRAUMATOLOGY | 26 |
| 7 | EXTENSION / SUPPORT ACCESSORIES | 33 |
| 8 | GENERAL ACCESSORIES | 43 |
| 9 | TRAYS | 48 |
| 10 | CARTS | 50 |
| 11 | NEURO | 52 |
| 12 | TABLE PADS | 53 |
| 13 | POWER CORDS | 59 |
| | Accessories List | 60 |

SIDE RAIL DIMENSIONS



actual size

NOTATION

MOT-VS600 Series : MOT-VS600D, MOT-VS600DH, MOT-VS600DK, MOT-VS600DHK

MOT-VS500 Series : MOT-VS500DK Series, MOT-VS500SK Series

MOT-VS500DK Series: MOT-VS500DK, MOT-VS500DK-IF

MOT-VS500SK Series : MOT-VS500SK, MOT-VS500SK-IF

MOT-5602 Series : MOT-5602SRMB, MOT-5602BW

MST-7201 Series : MST-7201B, MST-7201BX

MST-7300 Series : MST-7300B, MST-7300BX

R1 : Clamp for US Side Rail

R2 : Clamp for EU Side Rail

C1 : Clamp for US Side Rail. Size/Shape of accessory differs from the original.

C2 : Clamp for EU Side Rail. Size/Shape of accessory differs from the original.

NC: No Clamp. Requires additional clamp(s).

UE: No Clamp. Product size/shape differs from the original.

Requires additional clamp(s).


NOTE: The medical device registration status of each product on this catalog may be different according to each country's law. Please contact your Mizuho sales representative for additional information.

Non CE marked : Products with this mark are not confirmed their conformity to the provision of EU legislation MDR(EU) 2017/745. This mark does not apply to products which do not fall under medical devices legislation.


Locking Type Armboard



- Locks to side rail to prevent slipping
- Can be swiveled horizontally 180°



08-065-70 (mm)



08-065-80
08-065-85 (mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------|------|-----------|-----------------------|
| 08-065-70 | Locking Type Armboard | Each | JP | Pad not included |
| 08-065-80 | Locking Type Armboard (C1) | Each | US | |
| 08-065-85 | Locking Type Armboard (C2) | Each | EU | |

Pad


| Product Code | Description | Side Rail | Product Code | Description | Side Rail |
|--------------|-------------------------|-----------|--------------|---------------------------------|-----------|
| 08-065-08 | 40mm Pad for Armboard | JP | | | |
| 08-065-40 | 60mm Pad for Armboard | JP | 08-117-07 | 60mm Pad for Armboard (C1/C2) | US/EU |
| 08-065-45 | 100mm Pad for Armboard* | JP | 08-117-08 | 100mm Pad for Armboard (C1/C2)* | US/EU |

*Suitable when Radiographic Table Top is used with Pad

Double Axis Armboard



- Provides additional 15° or 30° vertical angulation at the elbow joint
- Prevents hyperextension of the arm
- Slits in armboard help to secure arm strap in place
- Locks to side rail to prevent slipping



(mm)




*Arm strap not included

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------|------|-----------|-----------------------|
| 08-065-51 | Double Axis Armboard | Each | JP | Pad not included |
| 08-065-51-R1 | Double Axis Armboard (R1) | Each | US | |
| 08-065-51-R2 | Double Axis Armboard (R2) | Each | EU | |

Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------------|------|-----------|---|
| 08-065-52 | 60mm Pad for Double Axis Armboard | Each | N/A | Suitable when Radiographic Table Top is used with pad |
| 08-065-53 | 100mm Pad for Double Axis Armboard | Each | N/A | |

Arm Strap




- Detachable
- Adjustable size
- Soft patient interface
- Washable

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------|------|-----------|-------------------------|
| 08-084-03 | Arm Strap | Each | N/A | L: 700 (mm), W: 80 (mm) |

IVR Armboard



- Armboard for IVR (Interventional Radiology)
- Radiolucent
- Clear imaging for catheter insertion
- Slides under table pad and is held in place by patient weight




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------|------|-----------|-----------------------|
| 08-065-72 | IVR Armboard | Each | N/A | Pad included |

Raised Armboard



- Adjusts to various patient body types and pad thicknesses
- For use in lateral and supine position
- Ball joint adjustability




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------------|------|-----------|---|
| 08-066-06 | Raised Armboard with Clamp (730mm) | Each | JP | |
| 08-066-10-UE | Raised Armboard (UE) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Armboard for Prone Position



- For use in prone position




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------------|------|-----------|---|
| 08-066-05-NC | Armboard for Prone Position (NC) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Lower Armboard



- For use in lateral, beach chair, and neuro park bench positions
- Ball joint adjustability
- Adjustable height
- Wide range of motion




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------|------|-----------|------------------------|
| 08-086-15 | Lower Armboard with Clamp | Each | JP | Pad and clamp included |
| 08-086-15-R1 | Lower Armboard (R1) | Each | US | |
| 08-086-15-R2 | Lower Armboard (R2) | Each | EU | |

Imaging Armboard



- Radiolucent top for C-arm access
- Legless design
- Ergonomic access to surgical site
- Securely attaches to any part of side rail



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------|------|-----------|-----------------------|
| 08-065-13 | Imaging Armboard | Each | JP | Pad not included |
| 08-065-13-R1 | Imaging Armboard (R1) | Each | US | |
| 08-065-13-R2 | Imaging Armboard (R2) | Each | EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------------|------|-----------|---|
| 08-065-26 | 50mm Pad for Imaging Armboard | Each | N/A | Suitable when Radiographic Table Top is used with pad |
| 08-065-27 | 90mm Pad for Imaging Armboard | Each | N/A | |

Arm Surgery Board



- Radiolucent top for C-arm access
- Securely attaches to any part of side rail
- Padless design




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|---------------------------------|
| 08-086-05 | Arm Surgery Board | Each | JP | Equipped with a side rail mount |
| 08-086-05-C1 | Arm Surgery Board with Side Rail Mount (C1) | Each | US | |
| 08-086-05-C2 | Arm Surgery Board with Side Rail Mount (C2) | Each | EU | |

Radiolucent Acrylic Armboard



- Slides under table pad and is held in place by patient weight and support rod
- Unique design allows for ergonomic access to the surgical site
- Added stability with adjustable support rod
- Folding design allows for easy storage
- Padless design



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------|------|-----------|-----------------------|
| 03-001-74 | Radiolucent Acrylic Armboard | Each | N/A | |

Arm Shield



- Slides under table pad and is held in place by patient weight
- Patient's arm is protected and secured in position



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information | W | D | H |
|--------------|-----------------------------|------|-----------|-----------------------|-----|-----|-----|
| 08-086-27 | Arm Shield with Pad (130mm) | Pair | N/A | Pad included | 400 | 250 | 130 |
| 08-086-28 | Arm Shield with Pad (160mm) | Pair | N/A | | 400 | 250 | 160 |
| 08-086-29 | Arm Shield with Pad (190mm) | Pair | N/A | | 400 | 250 | 190 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------|------|-----------|--|
| 08-069-45 | Pad for Arm Shield (150mm) | Pair | N/A | <ul style="list-style-type: none"> • Replacement pad • Pad thickness: 25mm |

Wrist Strap



- Holds patient's arm against the body
- Simple to use
- Used by bending over patient's forearm




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------|------|-----------|-----------------------|
| 08-079-00 | Wrist Strap (240mm) | Each | JP | |
| 08-079-00-R1 | Wrist Strap (240mm) (R1) | Each | US | |
| 08-079-00-R2 | Wrist Strap (240mm) (R2) | Each | EU | |

Shoulder Supports



- Prevents patient shifting when in Trendelenburg
- Easy to clean




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------|------|-----------|---|
| 08-068-04 | Shoulder Supports | Pair | JP,US,EU | <ul style="list-style-type: none"> • Clamp not included • Requires two Clamps for Flat Bar (Refer to p.45) JP : 08-111-11 US : 08-117-10 EU : 08-117-13 Universal : 08-111-20 |

Upper Arm Support II



- Adjusts patient's arm position freely
- Support base is radiolucent
- For use in supracondylar fracture surgery of the humerus
- Designed to keep the support base away from the post and avoid C-arm interference




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------|------|-----------|---|
| 08-064-22-NC | Upper Arm Support II (NC) | Each | JP,US,EU | <ul style="list-style-type: none"> Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Upper Arm Support Joint Type II



- Adjusts patient's arm position freely
- Adjustable in height and angle
- Support base is radiolucent
- For use in supracondylar fracture surgery of the humerus
- Designed to keep the support base away from the post and avoid C-arm interference



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------------------|------|-----------|---|
| 08-064-23-NC | Upper Arm Support Joint Type II (NC) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Levigator Stirrup II



- Available in a range of boot sizes to accommodate different patient types
- Integrated clamps fit instantly onto side rails of table
- Ergonomic grab-handle design and gas cylinder lift assist provides easier and smoother operation
- Light-weight design allows for quick and safe patient positioning
- Contoured boots with soft pads protect the entire lower leg to reduce the risk of superficial fibular nerve damage
- Lithotomy range : -35° low to 80° high/
10° adduction and 30° abduction



Levigator Stirrup II Set Main Unit + Boot + Pad

| | Product Code | Description | Unit | Side Rail | Technical Information |
|---|--------------|--|------|-----------|---|
| R | 18-070-91 | Levigator Stirrup II Set, Regular, Right Leg | Set | JP,US,EU | Component: 18-070-93 18-070-21 18-070-41 |
| L | 18-070-92 | Levigator Stirrup II Set, Regular, Left Leg | Set | JP,US,EU | Component: 18-070-94 18-070-22 18-070-42 |

*18-070-91 and 18-070-92 should be purchased together for both legs.


Main Unit



Main Unit

| | Product Code | Description | Unit | Side Rail |
|---|--------------|---|------|-----------|
| R | 18-070-93 | Levigator Stirrup II Main Unit, Right Leg | Each | JP,US,EU |
| L | 18-070-94 | Levigator Stirrup II Main Unit, Left Leg | Each | JP,US,EU |

*Right and left legs sold separately.




(mm)

Short Boot Extender

| | Product Code | Description | Unit | Side Rail |
|--|--------------|---------------------|------|-----------|
| | 18-070-50 | Short Boot Extender | Each | JP,US,EU |

*Order two pieces for both legs.

Non CE marked



- Ideal for patients whose feet do not reach the short boots.



Regular

Standard adult size



Boot

| | Product Code | Description | Unit |
|---|--------------|---|------|
| R | 18-070-21 | Boot for Levitator Stirrup II, Regular, Right Leg | Each |
| L | 18-070-22 | Boot for Levitator Stirrup II, Regular, Left Leg | Each |



Pad

| | Product Code | Description | Unit |
|---|--------------|---|------|
| R | 18-070-41 | Boot Pad for Levitator Stirrup II, Regular, Right Leg | Each |
| L | 18-070-42 | Boot Pad for Levitator Stirrup II, Regular, Left Leg | Each |

Fin Type

Lateral knee support further reduces risk of nerve damage



Boot

| | Product Code | Description | Unit |
|---|--------------|--|------|
| R | 18-070-23 | Boot for Levitator Stirrup II, Fin Type, Right Leg | Each |
| L | 18-070-24 | Boot for Levitator Stirrup II, Fin Type, Left Leg | Each |



Pad

| | Product Code | Description | Unit |
|---|--------------|--|------|
| R | 18-070-43 | Boot Pad for Levitator Stirrup II, Fin Type, Right Leg | Each |
| L | 18-070-44 | Boot Pad for Levitator Stirrup II, Fin Type, Left Leg | Each |

Wide Type

For large, overweight and obese patients



Boot

| | Product Code | Description | Unit |
|---|--------------|--|------|
| R | 18-070-27 | Boot for Levitator Stirrup II, Wide, Right Leg | Each |
| L | 18-070-28 | Boot for Levitator Stirrup II, Wide, Left Leg | Each |



Pad

| | Product Code | Description | Unit |
|---|--------------|--|------|
| R | 18-070-47 | Boot Pad for Levitator Stirrup II, Wide, Right Leg | Each |
| L | 18-070-48 | Boot Pad for Levitator Stirrup II, Wide, Left Leg | Each |

Short Type

For smaller patients



Boot

| | Product Code | Description | Unit |
|---|--------------|---|------|
| R | 18-070-25 | Boot for Levitator Stirrup II, Short, Right Leg | Each |
| L | 18-070-26 | Boot for Levitator Stirrup II, Short, Left Leg | Each |



Pad

| | Product Code | Description | Unit |
|---|--------------|---|------|
| R | 18-070-45 | Boot Pad for Levitator Stirrup II, Short, Right Leg | Each |
| L | 18-070-46 | Boot Pad for Levitator Stirrup II, Short, Left Leg | Each |


Levigator Cart



- Dedicated levigator stirrups cart with wheel castors
- For easy transport and storage



*Accessories in the image not included




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------|------|-----------|-----------------------|
| 18-070-61 | Levigator Stirrups Cart | Each | N/A | |

Transfer Board II



- For preoperative positioning of legs for lithotomy position
- Ultra-lightweight design
- Quick and easy mounting and removal




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 08-070-57 | Transfer Board II (Rail Fixation Type) | Each | JP,US,EU | |

Knee Crutches



- Hold patient's legs
- For use in lithotomy position



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|--|
| 08-070-00-UE | Knee Crutches (UE) | Pair | JP,US,EU | <ul style="list-style-type: none"> • Pad not included • Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|-----------------------|
| 08-070-13 | Urethane Pads for Knee Crutches | Pair | N/A | |

Pediatric Knee Crutches



- Hold pediatric patient's legs
- For use in lithotomy position



(mm)

Large

(mm)

| Product Code | Description | Unit | Side Rail | Technical Information | W | D |
|--------------|-------------------------------------|------|-----------|--|-----|-----|
| 08-070-16-NC | Pediatric Knee Crutches, Large (NC) | Pair | JP,US,EU | Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 | 100 | 150 |

Middle

(mm)

| Product Code | Description | Unit | Side Rail | Technical Information | W | D |
|--------------|--------------------------------------|------|-----------|--|----|-----|
| 08-070-17-NC | Pediatric Knee Crutches, Middle (NC) | Pair | JP,US,EU | Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 | 80 | 120 |

Small


(mm)

| Product Code | Description | Unit | Side Rail | Technical Information | W | D |
|--------------|-------------------------------------|------|-----------|--|----|----|
| 08-070-18-NC | Pediatric Knee Crutches, Small (NC) | Pair | JP,US,EU | Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 | 50 | 85 |

Knee Rests



- Provides a platform for kneeling
- Attached to leg section(s)
- Can be used with solid or split leg tables
- Can be used as leg extensions



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------------|------|-----------|---|
| 08-072-01 | Knee Rests | Pair | JP,US,EU | • Pads included • For the use with Radiographic Table Top, 08-072-02 will be applied • Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |
| 08-072-02 | Radiographic Tops for Knee Rests | Pair | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------|------|-----------|-----------------------|
| 08-072-03 | 70mm Pads for Knee Rests | Pair | N/A | Replacement pad |

Foot Rests



- Prevents patient from shifting when in reverse Trendelenburg
- Easy to clean




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|---|
| 08-071-03 | Foot Rests (320mm) | Pair | JP,US,EU | • Pad not included • Requires two Clamps for Flat Bar (Refer to p.45) JP : 08-111-11 US : 08-117-10 EU : 08-117-13 Universal : 08-111-20 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------------------|------|-----------|-----------------------|
| 08-071-05 | Mizuho Original Pads for Foot Rests | Pair | N/A | |

Head Rest



- Provides support of the patient's head in supine and prone positions



U-Shaped Head Rest

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------|------|-----------|-----------------------|
| 08-073-03 | U-Shaped Head Rest, L | Each | N/A | |
| 08-073-04 | U-Shaped Head Rest, S | Each | N/A | |




- Device for attaching U-shaped Head Rest / Mizuho Head Holder to the OR table
- Ball joint allows for easy head positioning


Table Attachment

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------|------|-----------|------------------------------------|
| 07-952-00 | Table Attachment | Each | N/A | Range of motion of pins: 290-390mm |

Narrow Head Section



- Improved access to patient head area



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------|------|-----------|---|
| 08-075-00 | Narrow Head Section | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MST-7201 Series MST-7300 Series MOT-3602 MOT-3502 MOT-6302 MOT-1602BW |
| 08-075-00-R1 | Narrow Head Section (R1) | Each | US | |
| 08-075-00-R2 | Narrow Head Section (R2) | Each | EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------------------|------|-----------|--|
| 08-083-69 | 70mm Pad for Narrow Head Section | Each | N/A | |
| 08-083-70 | 110mm Pad for Narrow Head Section | Each | N/A | *When used with Radiographic Table Top |

Head Supporter



- Provides secure head immobilization when in tilt and Trendelenburg
- Supports anesthetic management
- Multi-planar adjustment for optimal patient fit
- Integrated clamp attaches to side rail on either side of the table




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------|------|-----------|-----------------------|
| 08-075-06 | Head Supporter | Each | JP | |
| 08-075-06-R1 | Head Supporter (R1) | Each | US | |
| 08-075-06-R2 | Head Supporter (R2) | Each | EU | |

Adapter for Carbon Fiber Extension



- To be used to attach head frame to carbon fiber board



(mm)




*Table attachment and head holder not included

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------------|------|-----------|---|
| 08-075-15 | Adapter for Carbon Fiber Extension | Each | JP | Applicable table and accessories: VACB-3006B Carbon Fiber Extension II (08-070-54) (Refer to p.39) |

Attachment for Mayfield Head Support



- To be used to attach base unit of skull clamp etc. to table




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------------------|------|-----------|--|
| 08-075-11 | Attachment for Mayfield Head Support | Each | JP,US,EU | Applicable tables: MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MOS-1300B/1302B MST-7201Series MST-7300 Series SPL-331 MOT-3602 MOT-3502 MOT-6302 |

Chest & Waist Support



- May be used in supine or lateral to support body
- Easy positioning with adjustable height and ball joint




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------|------|-----------|--|
| 08-069-15-UE | Chest & Waist Support (UE) | Each | JP,US,EU | Requires Clamp for Flat Bar (Refer to p.49) JP : 08-111-11 US : 08-117-10 EU : 08-117-13 Universal : 08-111-20 |

Lateral Support



- Provides a wide range of adjustment so that operative field is not obstructed
- Secure support in lateral and other positions
- Easy positioning with adjustable height and ball joint




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------|------|-----------|--|
| 08-069-24-NC | Lateral Support (NC) | Each | JP,US,EU | Requires Clamp for Flat Bar (Refer to p.45) JP : 08-111-11 US : 08-117-10 EU : 08-117-13 Universal : 08-111-20 |

Adjustable Body Support



- Hydraulic locking mechanism
- Wide range of adjustment
- Holds patient rigidly even in steep lateral tilt
- Easy positioning with adjustable height and ball joint



(mm)



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------|------|-----------|---|
| 08-080-11-NC | Adjustable Body Support (NC) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Body Strap




- Detachable
- Adjustable size
- Soft patient interface



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------|------|-----------|---------------------------------------|
| 08-084-02 | Body Strap (1700mm) | Each | JP,US,EU | L: 1700 (mm), W: 100 (mm) Washable |

Body Strap with Support Pad

For quickly and securely fixing the body



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------|------|-----------|-----------------------|
| 08-084-00 | Leg Strap | Each | JP,US,EU | L: 1890 (mm) |
| 08-084-01 | Body Strap (2500mm) | Each | JP,US,EU | L: 2500 (mm) |

Radiographic Table Top / Fixed Type




- Designed for optimal X-ray cassette use

Choose required sections or all sections for whole table configuration from item list (pp.20-21).


Note: Table tops for MOT-5602BW and MST-7300 Series are sold as a set respectively

MOT-VS700UKI with Double Leg Section (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-092-20 | Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS700UKI with Single Leg Section (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-092-20 | Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-090-60 | Radiographic Table Tops Fixed Type for Single Leg Section Type A | Each |

MOT-VS600D / DH (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-090-10 | Radiographic Table Tops Fixed Type for Back Section Type A | Each |
| 3 18-090-20 | Radiographic Table Tops Fixed Type for Seat Section Type A | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS600DK / DHK (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-092-20 | Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS500DK (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-092-21 | Radiographic Table Tops Fixed Type for Seat Section Type KB | Each |
| 4 18-090-51 | Radiographic Table Tops Fixed Type for Double Leg Section Type B | Pair |

MOT-VS500SK (Fixed Type)




| Product Code | Description | Unit |
|--------------|---|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-092-21 | Radiographic Table Tops Fixed Type for Seat Section Type KB | Each |
| 4 18-090-61 | Radiographic Table Tops Fixed Type for Single Leg Section Type A-II | Each |

MOT-VE100D (Fixed Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-090-10 | Radiographic Table Tops Fixed Type for Back Section Type A | Each |
| 3 18-090-20 | Radiographic Table Tops Fixed Type for Seat Section Type A | Each |
| 4 18-090-51 | Radiographic Table Tops Fixed Type for Double Leg Section Type B | Pair |

MOT-5602BW (Fixed Type)



| Product Code | Description | Unit |
|--------------|---|------|
| 08-089-72 | Radiographic Table Tops Fixed Type for 5602BW | Set |

MST-7300 Series (Fixed Type)



| Product Code | Description | Unit |
|--------------|--|------|
| 08-089-76 | Radiographic Table Tops Fixed Type for 7300 Series | Set |

Radiographic Table Top / Fixed Type for Digital Cassette




- Designed for optimal X-ray cassette use
- The digital radiographic table top fixed type has a wide width and the digital cassette can be inserted easily

Choose required sections or all sections for whole table configuration from item list (pp.22-23).


Note: Table tops for MOT-5602BW and MST-7300 Series are sold as a set respectively

MOT-VS700UKI with Double Leg Section (Digital Type)




| Product Code | Description | Unit |
|--------------|---|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-093-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS700UKI with Single Leg Section (Digital Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-093-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-091-60 | Digital Radiographic Table Tops Fixed Type for Single Leg Section Type A | Each |

MOT-VS600D / DH (Digital Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-091-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type A | Each |
| 3 18-091-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type A | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS600DK / DHK (Digital Type)




| Product Code | Description | Unit |
|--------------|---|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-093-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KA | Each |
| 4 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | Pair |

MOT-VS500DK (Digital Type)




| Product Code | Description | Unit |
|--------------|---|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-093-21 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KB | Each |
| 4 18-090-51 | Radiographic Table Tops Fixed Type for Double Leg Section Type B | Pair |

MOT-VS500SK (Digital Type)




| Product Code | Description | Unit |
|--------------|---|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | Each |
| 3 18-093-21 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KB | Each |
| 4 18-091-61 | Digital Radiographic Table Tops Fixed Type for Single Leg Section Type A-II | Each |

MOT-VE100D (Digital Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 1 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | Each |
| 2 18-091-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type A | Each |
| 3 18-091-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type A | Each |
| 4 18-090-51 | Radiographic Table Tops Fixed Type for Double Leg Section Type B | Pair |

MOT-5602BW (Digital Type)



| Product Code | Description | Unit |
|--------------|---|------|
| 08-089-82 | Digital Radiographic Table Tops Fixed Type for 5602BW | Set |

MST-7300 Series (Digital Type)




| Product Code | Description | Unit |
|--------------|--|------|
| 08-089-86 | Digital Radiographic Table Tops Fixed Type for 7300 Series | Set |

Radiographic Table Top / Removable Type




- Designed for optimal X-ray cassette use
- Refer to the list below for applicable type for your table
- Below models are sold as a set

MOT-5602 Series (Removable type)




| Product Code | Description | Unit |
|--------------|---|------|
| 08-088-20 | Radiographic Table Top for MOT-5602SRMB | Set |

MOS-1302B (Removable type)




| Product Code | Description | Unit |
|--------------|-----------------------------------|------|
| 08-088-59 | Radiographic Table Tops for 1302B | Set |

SPL-331 (Removable type)




| Product Code | Description | Unit |
|--------------|-------------------------------------|------|
| 08-088-22 | Radiographic Table Tops for SPL-331 | Set |

Cassette Tray



- Holds x-ray cassettes with a hook-type handle




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------|------|-----------|-----------------------|
| 08-094-00 | Cassette Tray (390 x 466 mm) | Each | N/A | |

Adjustable Cassette Tray



- Various cassettes can be used by adjusting width of the tray
- Scaled handle indicates insertion position
- Handle is removable



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------|------|-----------|-----------------------|
| 08-095-50 | Adjustable Cassette Tray | Each | N/A | |



Leg Traction Device MOC-142



- Free standing lower extremity traction device
- Rolled into position and locked in place
- Does not attach to table
- Can be used with any MOT table including rotating top table

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------------|------|-----------|-----------------------|
| 08-051-00 | Leg Traction Device MOC-142 | Set | N/A | |

Non CE marked



Traction Device

- Applicable surgery: femoral fracture, tibia fracture
- Good C-arm access

| Applicable Tables | Technical Information |
|---|--|
| MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MOT-3602 MOT-3502 | • Requires V Attachment and Cart for Traction Arm compatible with operating table (Refer to p.29 and p.30) |

Recommended Combination of Accessories

Femoral Surgery / Supine



| No. | Product Code | Description | Unit | Side Rail |
|-----|--------------|---|------|-----------|
| 1 | 08-115-01 | Traction Device Arm | Each | N/A |
| 2 | 08-115-02 | Extension Arm, Right/Left | Pair | N/A |
| 3 | 08-115-04 | Traction Device Right/Left | Pair | N/A |
| 4 | 08-115-05 | Traction Shoe Right/Left | Pair | N/A |
| 5 | 08-115-07 | Inner Thigh Column with Pad | Each | N/A |
| 6 | 08-115-27 | Support Rods | Pair | N/A |
| 7 | 08-115-76 | Wide Pelvic Support II | Each | N/A |
| 8 | 08-115-77 | Transfer Board II (Pin Fixation Type) ^{*1} | Each | N/A |
| 9 | 08-115-10 | Lower Limb Support ^{*1} | Each | N/A |

*1 Not shown in picture.

Femoral Surgery / Lateral



| No. | Product Code | Description | Unit | Side Rail |
|-----|--------------|---|------|-----------|
| 1 | 08-115-01 | Traction Device Arm | Each | N/A |
| 2 | 08-115-02 | Extension Arm, Right/Left | Pair | N/A |
| 3 | 08-115-04 | Traction Device Right/Left | Pair | N/A |
| 4 | 08-115-05 | Traction Shoe Right/Left | Pair | N/A |
| 6 | 08-115-27 | Support Rods | Pair | N/A |
| 9 | 08-115-10 | Lower Limb Support | Each | N/A |
| 10 | 08-115-11 | Inner Thigh Support with Pad for Traction in Lateral Posture | Each | N/A |
| 11 | 08-115-12 | Lower Thigh Support for Traction in Lateral Posture ^{*1} | Each | N/A |
| 12 | 08-115-26 | Pelvic Support Board ^{*2} | Each | JP |
| | 08-115-26-R1 | Pelvic Support Board (R1) ^{*2} | Each | US |
| | 08-115-26-R2 | Pelvic Support Board (R2) ^{*2} | Each | EU |

*1 Not shown in picture. *2 Please select to match the side rail for your table.

Tibia Surgery



| No. | Product Code | Description | Unit | Side Rail |
|-----|--------------|---|------|-----------|
| 1 | 08-115-01 | Traction Device Arm | Each | N/A |
| 2 | 08-115-02 | Extension Arm, Right/Left | Pair | N/A |
| 3 | 08-115-04 | Traction Device Right/Left | Pair | N/A |
| 4 | 08-115-05 | Traction Shoe Right/Left | Pair | N/A |
| 6 | 08-115-27 | Support Rods | Pair | N/A |
| 9 | 08-115-10 | Lower Limb Support | Each | N/A |
| 12 | 08-115-26 | Pelvic Support Board ^{*2} | Each | JP |
| | 08-115-26-R1 | Pelvic Support Board (R1) ^{*2} | Each | US |
| | 08-115-26-R2 | Pelvic Support Board (R2) ^{*2} | Each | EU |
| 13 | 08-115-09 | L-Shaped Knee Rest with Pad | Each | N/A |

*2 Please select to match the side rail for your table.

Optional Accessories


| No. | Product Code | Description | Unit | Side Rail |
|-----|--------------|---|------|-----------|
| 14 | 08-081-07 | Accessory Cart for Traction Device | Each | JP |
| | 08-081-07-R1 | Accessory Cart for Traction Device (R1) | Each | US |
| | 08-081-07-R2 | Accessory Cart for Traction Device (R2) | Each | EU |
| 15 | 08-115-08 | Pelvic Support for Pediatric | Each | N/A |
| 16 | 08-115-17 | Traction Pulley | Each | N/A |
| 17 | 08-115-18 | Direct Traction Device | Each | N/A |
| 18 | 08-115-19 | Direct Traction Hook Attachment | Each | N/A |
| 19 | 08-115-20 | Movable Pelvic Support | Each | N/A |
| 20 | 08-115-36 | Inner Thigh Column | Each | N/A |
| 21 | 08-115-37 | Pad for Inner Thigh Column | Each | N/A |
| 22 | 08-115-38 | Lower Limb Support Clamp | Each | N/A |
| 23 | 08-115-43 | Carbon Fiber L Shaped Knee Rest with Pad | Each | N/A |
| 24 | 08-115-45 | Direct Traction Device for Borshards | Each | N/A |
| 25 | 08-115-73 | Lower Limb Support without Fixing Bracket | Each | N/A |

Traction Accessories



- 1 Traction Device Arm**
Product Code 08-115-01
- 2 Extension Arm, Right/Left**
Product Code 08-115-02
- 3 Traction Device Right/Left**
Product Code 08-115-04
- 4 Traction Shoe Right/Left**
Product Code 08-115-05

- 5 Inner Thigh Column with Pad**
Product Code 08-115-07
- 6 Support Rods**
Product Code 08-115-27
- 7 Wide Pelvic Support II**
Product Code 08-115-76
- 8 Transfer Board II (Pin Fixation Type)**
Product Code 08-115-77



9 Lower Limb Support
Product Code 08-115-10



12 Pelvic Support Board
Product Code 08-115-26



10 Inner Thigh Support with Pad for Traction in Lateral Posture
Product Code 08-115-11

12 Pelvic Support Board (R1)
Product Code 08-115-26-R1

12 Pelvic Support Board (R2)
Product Code 08-115-26-R2



11 Lower Thigh Support for Traction in Lateral Posture
Product Code 08-115-12



13 L-Shaped Knee Rest with Pad
Product Code 08-115-09

Optional Accessories



14 Accessory Cart for Traction Device
Product Code 08-081-07



16 Traction Pulley
Product Code 08-115-17



19 Movable Pelvic Support
Product Code 08-115-20



23 Carbon Fiber L Shaped Knee Rest with Pad
Product Code 08-115-43

14 Accessory Cart for Traction Device (R1)
Product Code 08-081-07-R1



17 Direct Traction Device
Product Code 08-115-18



20 Inner Thigh Column
Product Code 08-115-36



24 Direct Traction Device for Borshards
Product Code 08-115-45

14 Accessory Cart for Traction Device (R2)
Product Code 08-081-07-R2




15 Pelvic Support for Pediatric
Product Code 08-115-08



18 Direct Traction Hook Attachment
Product Code 08-115-19

21 Pad for Inner Thigh Column
Product Code 08-115-37



22 Lower Limb Support without Fixing Bracket
Product Code 08-115-38




25 Lower Limb Support without Fixing Bracket
Product Code 08-115-73

Cart for Traction Arm



- Simplifies installation of Leg Traction Device




| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|--------------------------------|------|-----------|---|
| 08-115-21 | Cart for Traction Arm | Each | N/A | MOS-1302B |
| 08-115-39 | Cart for Traction Arm Type II | Each | N/A | MOT-5602 Series |
| 08-115-40 | Cart for Traction Arm Type III | Each | N/A | MOT-5602 Series *When used with Radiographic Table Tops |
| 08-115-41 | Cart for Traction Arm Type IV | Each | N/A | MOT-3502 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MOT-3602 *When used with Radiographic Table Tops |
| 08-115-42 | Cart for Traction Arm Type V | Each | N/A | MOT-3502 *When used with Radiographic Table Tops |
| 08-115-55 | Cart for Traction Arm Type VI | Each | N/A | MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MOT-3602 |

V Attachment II Side Rail Type



- Required to attach the Traction Device to the operating table
- Built-in safety hook for increased security and stability
- Attaches to side rail of seat section




| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|--|------|-----------|---|
| 08-115-80 | V Attachment II Side Rail Type | Each | JP | MOT-5602 Series MOT-VS700UKI |
| 08-115-80-R1 | V Attachment II Side Rail Type (R1) | Each | US | MOT-VS600 Series |
| 08-115-80-R2 | V Attachment II Side Rail Type (R2) | Each | EU | MOT-VS500DK Series MOT-VE100D |
| 08-115-81 | V Attachment II Side Rail Type for Radiographic Table Top | Each | JP | MOT-5602 Series MOT-VS700UKI |
| 08-115-81-R1 | V Attachment II Side Rail Type Radiographic Table Top (R1) | Each | US | MOT-VS600 Series MOT-VS500DK Series |
| 08-115-81-R2 | V Attachment II Side Rail Type Radiographic Table Top (R2) | Each | EU | MOT-VE100D *When used with Radiographic Table Tops |

V Attachment Pin Type



- Required to attach the Traction Device to the operating table
- After removing the leg section, attach it by inserting the pin of the seat section into the attachment hole




| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|---|------|-----------|---|
| 08-115-23 | V Attachment Pin Type | Each | JP | MOT-VS500SK Series |
| 08-115-23-R1 | V Attachment Pin Type (R1) | Each | US | MOT-3502 |
| 08-115-23-R2 | V Attachment Pin Type (R2) | Each | EU | MOT-3602 |
| 08-115-24 | V Attachment Pin Type for Radiographic Table Top | Each | JP | MOT-VS500SK Series |
| 08-115-24-R1 | V Attachment Pin Type for Radiographic Table Top (R1) | Each | US | MOT-3502 |
| 08-115-24-R2 | V Attachment Pin Type for Radiographic Table Top (R2) | Each | EU | MOT-3602 *When used with Radiographic Table Tops |

Spinal Frame



- For spinal surgery
- Curved frame can be flexed with a handle
- Enables space between vertebrae




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|-----------------------|
| 02-120-00 | GM-98 Spinal Frame | Each | N/A | |

Non CE marked


Four-Point Spinal Frame



- For spinal surgery
- Supports patient at four points with soft polymer gel cushions
- Reduced abdominal pressure by allowing patient's abdomen to hang freely




(mm)




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------|------|-----------|-----------------------|
| 02-121-00 | Four-Point Spinal Frame | Each | N/A | Pillow included |


Radiolucent Carbon Fiber Four-Point Spinal Frame



- For spinal surgery
- 360° radiolucent
- Supports patient at four points with soft polymer gel cushions
- Reduced abdominal pressure by allowing patient's abdomen to hang freely



(mm)




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 02-121-10 | Radiolucent Carbon Fiber Four Point Spinal Frame | Each | N/A | Pillow included |

Shoulder Positioner



- For arthroscopic, open shoulder surgery and joint replacement
- Provides secure beach chair position
- Can be lifted by electric power
- Back plate is detachable so that surgical site will be radiolucent
- Good C-arm access



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|--|
| 08-068-22 | Shoulder Positioner without Pad | Each | JP,US,EU | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series MST-7300 Series MOT-3602 |

Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------------|------|-----------|-----------------------|
| 08-068-23 | Pad for Shoulder Positioner | Set | N/A | |


Replacement Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 08-068-24 | Shoulder Positioner Replacement Strap for Head Support | Each | N/A | |
| 08-068-25 | Shoulder Positioner Replacement Pad, Head Section | Each | N/A | |
| 08-068-26 | Shoulder Positioner Replacement Pad, Shoulder Section | Each | N/A | |
| 08-068-27 | Shoulder Positioner Replacement Pad, Back Section | Each | N/A | |

Pubis Support



- Supports pubis in lateral procedure
- Rotatable pad enables more stable support by conforming to patient's body




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|--|
| 08-078-32-NC | Pubis Support (NC) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Table Side Extension



- Extends width of table by attaching to side rail



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------|------|-----------|-----------------------|
| 08-086-06 | Table Side Extension | Each | JP | Pad not included |
| 08-086-06-R1 | Table Side Extension (R1) | Each | US | |
| 08-086-06-R2 | Table Side Extension (R2) | Each | EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------------|------|-----------|---|
| 08-086-33 | 70mm Pad for Table Side Extension | Each | N/A | Suitable when Radiographic Table Top is used with pad |
| 08-086-35 | 110mm Pad for Table Side Extension | Each | N/A | |

Head Extension II



- Extends length of table to the head side
- Install between the head and back section
- Provides wider imaging area



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------|------|-----------|---|
| 08-075-90 | Head Extension II | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MST-7201 Series MST-7300 Series MOT-3602 MOT-3502 MOT-6302 MOT-1602BW |
| 08-075-90-R1 | Head Extension II (R1) | Each | US | |
| 08-075-90-R2 | Head Extension II (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|-----------------------|
| 08-083-02 | Radiographic Table Tops Fixed Type for Head Extension (240mm) | Each | JP,US,EU | |
| 08-083-22 | Digital Radiographic Table Tops Fixed Type for Head Extension (240mm) | Each | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|---|------|-----------|------------------------------|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-01 | SUSANO Pad for Head Section Type A | Each | N/A | Flame-retardant / Waterproof |
| Standard | | | | |
| 18-042-00 | Mizuho Original Pad for Head Section Type A | Each | N/A | |

Seat Extension




- Extends seat section by attaching to side rail
- Provides wider imaging area



08-075-36

(mm)



08-075-42

(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------|------|-----------|---|
| 08-075-36 | Seat Extension (350mm) | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5602 Series MOT-VS500DK Series MOT-VE100D MOT-1602BW • Transfer Board II (08-070-57) can be attached for holding patient's legs (Refer to p.12) |
| 08-075-36-R1 | Seat Extension (350mm) (R1) | Each | US | |
| 08-075-36-R2 | Seat Extension (350mm) (R2) | Each | EU | |
| 08-075-42 | Seat Extension for 5701 | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5701 MOT-VS600 Series • Leg section can be attached |
| 08-075-42-R1 | Seat Extension for 5701 (R1) | Each | US | |
| 08-075-42-R2 | Seat Extension for 5701 (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|--|
| 08-083-01 | Radiographic Table Tops Fixed Type for Seat Extension (350mm) | Each | JP,US,EU | To be used for: 08-075-36 08-075-36-R1 08-075-36-R2 |
| 08-083-21 | Digital Radiographic Table Tops Fixed Type for Seat Extension (350mm) | Each | JP,US,EU | |
| 08-083-00 | Radiographic Table Tops Fixed Type for Seat Extension (300mm) | Each | JP,US,EU | |
| 08-083-20 | Digital Radiographic Table Tops Fixed Type for Seat Extension (300mm) | Each | JP,US,EU | To be used for: 08-075-42 08-075-42-R1 08-075-42-R2 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|--|------|-----------|---|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-83 | SUSANO Pad for Seat Extension (350mm) | Each | N/A | <ul style="list-style-type: none"> • Flame-retardant / Waterproof • To be used for: 08-075-36 08-075-36-R1 08-075-36-R2 |
| 18-042-82 | SUSANO Pad for Seat Extension (300mm) | Each | N/A | <ul style="list-style-type: none"> • Flame-retardant / Waterproof • To be used for: 08-075-42 08-075-42-R1 08-075-42-R2 |
| Standard | | | | |
| 08-083-61 | Mizuho Original Pad for Seat Extension (350mm) | Each | N/A | To be used for: 08-075-36 08-075-36-R1 08-075-36-R2 |
| 08-083-60 | Mizuho Original Pad for Seat Extension (300mm) | Each | N/A | To be used for: 08-075-42 08-075-42-R1 08-075-42-R2 |

Seat Extension II



- Provides wider imaging area by extending the seat section



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------------|------|-----------|--|
| 08-075-52 | Seat Extension II (300mm) | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-VS700UKI • Leg section can be attached |
| 08-075-52-R1 | Seat Extension II (300mm) (R1) | Each | US | |
| 08-075-52-R2 | Seat Extension II (300mm) (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|--|
| 08-083-00 | Radiographic Table Tops Fixed Type for Seat Extension (300mm) | Each | JP,US,EU | To be used for: 08-075-52 08-075-52-R1 08-075-52-R2 |
| 08-083-20 | Digital Radiographic Table Tops Fixed Type for Seat Extension (300mm) | Each | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|--|------|-----------|---|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-82 | SUSANO Pad for Seat Extension (300mm) | Each | N/A | <ul style="list-style-type: none"> • Flame-retardant / Waterproof • To be used for: 08-075-52 08-075-52-R1 08-075-52-R2 |
| Standard | | | | |
| 08-083-60 | Mizuho Original Pad for Seat Extension (300mm) | Each | N/A | To be used for: 08-075-52 08-075-52-R1 08-075-52-R2 |

Leg Extension



- Extends leg section by attaching to side rail



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------------|------|-----------|--|
| 08-075-30 | Leg Extension for 5700 Series | Pair | JP | • Pad not included • Applicable tables: MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MST-7201 Series MST-7300 Series MOT-3602 MOT-3502 MOT-6302 MOT-1602BW SPL-331 |
| 08-075-30-R1 | Leg Extension for 5700 Series (R1) | Pair | US | |
| 08-075-30-R2 | Leg Extension for 5700 Series (R2) | Pair | EU | |


Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 08-083-04 | Radiographic Table Tops Fixed Type for Leg Extension | Pair | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|---------------------------------------|------|-----------|------------------------------|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-84 | SUSANO Pad for Leg Extension | Pair | N/A | Flame-retardant / Waterproof |
| Standard | | | | |
| 08-083-62 | Mizuho Original Pad for Leg Extension | Pair | N/A | |

Side Rail Extension



- For attachment of additional accessories
- Easily attaches to side rails





| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------------------|------|-----------|-----------------------|
| 08-115-32 | Side Rail Extension, 220mm (Pair) | Pair | JP | |



Uro Drain Tray Set





- Collects waste fluids during urology and gynecology procedures
- Stainless steel tray stowed under table when not in use




Funnel Type A





Funnel Type B



Funnel Type C



Funnel Type E



Waste Container,
Max 9000cc



Hose

(mm)


| Product Code | Description | Unit | Side Rail | Applicable Tables | Funnel |
|--------------|---|------|-----------|--|--------|
| 08-095-08 | Uro Drain Tray Set for 5602SRMB | Set | JP,US,EU | MOT-5602 Series | Type A |
| 08-095-13 | Uro Drain Tray Set for 3600 & 3500 Series, 5801 | Set | JP,US,EU | MOT-VS500 Series MOT-VE100D MOT-3602 MOT-3502 | Type B |
| 08-095-15 | Uro Drain Tray Set for 5701 | Set | JP,US,EU | MOT-5701 MOT-VS600 Series | Type C |
| 08-095-19 | Uro Drain Tray for SPL-331 | Set | JP | SPL-331 | Type E |
| 08-095-20 | Waste Container with Hose | Set | N/A | | |

*For SPL-331, the complete uro drain tray set will be the combination of 08-095-19 and 08-095-20.

Back Section Board Attachment for Pediatric



- 40cm width allows for improved surgeon ergonomics



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|--|
| 08-068-31 | Back Section Board Attachment for Pediatric | Each | JP | <ul style="list-style-type: none"> • Pad not included. • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500SK Series MST-7300 Series MOT-3602 *MOT-VS500DK Series *MOT-VE100D *When used with Leg-Side Reversal Adapter |
| 08-068-31-R1 | Back Section Board Attachment for Pediatric (R1) | Each | US | |
| 08-068-31-R2 | Back Section Board Attachment for Pediatric (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 08-083-06 | Radiographic Table Tops Fixed Type for Back Section Board Attachment for Pediatric | Each | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|---|------|-----------|------------------------------|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-81 | SUSANO Pad for Back Section Board Attachment for Pediatric | Each | N/A | Flame-retardant / Waterproof |
| Standard | | | | |
| 08-083-63 | Mizuho Original Pad for Back Section Board Attachment for Pediatric | Each | N/A | |

Back Section Board Short



- To be used to adjust table length according to patient size
- Shorter table length allows for easier anesthesia control



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------------|------|-----------|---|
| 08-007-21 | Back Section Board Short | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series |
| 08-007-21-R1 | Back Section Board Short (R1) | Each | US | |
| 08-007-21-R2 | Back Section Board Short (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|-----------------------|
| 08-083-03 | Radiographic Table Tops Fixed Type for Back Section Board Short | Each | JP,US,EU | |
| 08-083-23 | Digital Radiographic Table Tops Fixed Type for Back Section Board Short | Each | JP,US,EU | |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--|---|------|-----------|------------------------------|
| SUSANO - Flame-retardant / Waterproof | | | | |
| 18-042-80 | SUSANO Pad for Back Section Board Short | Each | N/A | Flame-retardant / Waterproof |
| Standard | | | | |
| 08-083-64 | Mizuho Original Pad (252mm) | Each | N/A | |

Kidney Bridge



- Enhances the flex of patient
- Can be lifted up to 270mm (including pad)




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------|------|-----------|---|
| 08-096-02 | Kidney Bridge with Pad | Each | JP | <ul style="list-style-type: none"> • Manual operation • Applicable tables: MOT-5602 Series MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500 Series MOT-VE100D MST-7201 Series MST-7300 Series SPL-331 MOT-3602 MOT-3502 MOT-6302 MOT-1602BW |

100cm Carbon Fiber Extension II



- Radiolucent
- 360-degree imaging
- Quick and easy attachment
- Length can be chosen according to patient size, type of surgery



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|---|
| 08-070-54 | 100cm Carbon Fiber Extension II | Each | JP,US,EU | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500SK Series MST-7300 Series MOT-3602 MOT-3502 MOT-6302 *MOT-VS500DK Series *MOT-VE100D *When used with Leg-Side Reversal Adapter • For Side Rail Adaptor refer to p.46 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|---|
| 08-083-65 | 70mm Mizuho Original Pad for 100cm Carbon Fiber Extension | Each | N/A | |
| 08-083-66 | 110mm Mizuho Original Pad for 100cm Carbon Fiber Extension | Each | N/A | Suitable when Radiographic Table Top is used with pad |

122cm Carbon Fiber Extension II



- Radiolucent
- 360-degree imaging
- Quick and easy attachment
- Length can be chosen according to patient size, type of surgery
- Added stability with adjustable support rod



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|---|
| 08-070-55 | 122cm Carbon Fiber Extension II | Each | JP,US,EU | <ul style="list-style-type: none"> • Pad not included • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500SK Series MST-7300 Series MOT-3602 MOT-3502 MOT-6302 *MOT-VS500DK Series *MOT-VE100D *When used with Leg-Side Reversal Adapter • For Side Rail Adaptor refer to p.46 |


Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|---|
| 08-083-67 | 70mm Mizuho Original Pad for 122cm Carbon Fiber Extension | Each | N/A | |
| 08-083-68 | 110mm Mizuho Original Pad for 122cm Carbon Fiber Extension | Each | N/A | Suitable when Radiographic Table Top is used with pad |

Leg-Side Reversal Adapter


- An adapter for reverse position
- After removing the leg section and attaching the adapter to the seat section, you can attach adaptive accessories

Leg-Side Reversal Adapter A




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-----------------------------|------|-----------|--|
| 08-077-70 | Leg-Side Reversal Adapter A | Each | JP,US,EU | <ul style="list-style-type: none"> • Pad not included • Applicable Table: MOT-VS700UKI |

Leg-Side Reversal Adapter B



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------------|------|-----------|--|
| 08-077-60 | Leg-Side Reversal Adapter B | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable Table: MOT-VS600 Series |
| 08-077-60-R1 | Leg-Side Reversal Adapter B (R1) | Each | US | |
| 08-077-60-R2 | Leg-Side Reversal Adapter B (R2) | Each | EU | |

Leg-Side Reversal Adapter C



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------------|------|-----------|---|
| 08-077-50 | Leg-Side Reversal Adapter C | Each | JP | <ul style="list-style-type: none"> • Pad not included • Applicable Tables: MOT-VS500DK Series, MOT-VE100D |
| 08-077-50-R1 | Leg-Side Reversal Adapter C (R1) | Each | US | |
| 08-077-50-R2 | Leg-Side Reversal Adapter C (R2) | Each | EU | |

Radiographic Table Top

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|---|
| 18-097-00 | Digital Radiographic Table Top for Leg-Side Reversal Adapter | Each | JP,US,EU | <ul style="list-style-type: none"> • To be used for: Leg-Side Reversal Adapter A, B, C |

Pad

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|---|
| 18-042-90 | Mizuho Original Pad for Leg-Side Reversal Adapter | Each | N/A | <ul style="list-style-type: none"> • To be used for: Leg-Side Reversal Adapter A, B, C |

Foot Switch

- For table operation
- Pedal function:
 - Left: Elevation
 - Middle: Trendelenburg
 - Right: Lateral tilt



| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|-------------|------|-----------|---------------------------------|
| 18-443-43 | Foot Switch | Each | N/A | MOT-VS700UKI MST-7300 Series |



| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|-------------------|------|-----------|--|
| 18-443-50 | Foot Switch TypeA | Each | N/A | MOT-VS600 Series MOT-VS500DK-IF MOT-VS500SK-IF |




| Product Code | Description | Unit | Side Rail | Applicable Tables |
|--------------|-------------------|------|-----------|--|
| 18-443-20 | Foot Switch TypeB | Each | N/A | MOT-5602 Series MOT-3502 MOT-6302 MOS-1300B/1302B VACB-3006B |

Anesthesia Screen L Type



- For use with sterile covers
- Adjustable height
- Adjustable horizontal orientation




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|-----------------------|
| 08-062-00 | Anesthesia Screen L Type with Clamp | Each | JP | |
| 08-062-00-R1 | Anesthesia Screen L Type with Clamp (R1) | Each | US | |
| 08-062-00-R2 | Anesthesia Screen L Type with Clamp (R2) | Each | EU | |

Anesthesia Screen Extension



- For use with sterile covers
- Adjustable height
- Adjustable width




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|-----------------------|
| 08-061-40 | Anesthesia Screen Extension with Clamp | Each | JP | |
| 08-061-40-R1 | Anesthesia Screen Extension with Clamp (R1) | Each | US | |
| 08-061-40-R2 | Anesthesia Screen Extension with Clamp (R2) | Each | EU | |

Anesthesia Screen Wing Type



- For use with sterile covers
- Adjustable height
- Widens support surface for sterile cover
- Wings can be swivelled



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------------------|------|-----------|---|
| 08-061-20-UE | Anesthesia Screen Wing Type (UE) | Each | JP,US,EU | Requires two Clamps for Flat Bar (Refer to p.45) JP : 08-111-11 US : 08-117-10 EU : 08-117-13 Universal : 08-111-20 |

Clamp



Clamp for Round Bar

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------------|------|-----------|-------------------------|
| 08-110-01 | Clamp for Round Bar | Each | JP | φ16mm |
| 08-117-09 | Clamp for Round Bar (R1) | Each | US | φ16mm |
| 08-117-11 | Clamp for Round Bar (R2) | Each | EU | φ16mm |
| 09-228-12 | Side Rail Clamp (φ18mm) | Each | JP | φ18mm Non C € marked |



Universal Clamp for Round Bar

Can be used for JP, US, EU side rails

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------------|------|-----------|-----------------------|
| 08-110-20 | Universal Clamp for Round Bar | Each | JP,US,EU | φ16mm |



Clamp for Flat Bar

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------------|------|-----------|-----------------------|
| 08-111-11 | Clamp for Flat Bar Ver.2 | Each | JP | |
| 08-117-10 | Clamp for Flat Bar Ver.2 (R1) | Each | US | |
| 08-117-13 | Clamp for Flat Bar Ver.2 (R2) | Each | EU | |



Universal Clamp for Flat Bar

Can be used for JP, US, EU side rails

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|------------------------------|------|-----------|-----------------------|
| 08-111-20 | Universal Clamp for Flat Bar | Each | JP,US,EU | |




Multi Clamp

Can be used for accessories with both flat and round bar

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------|------|-----------|-----------------------|
| 08-110-11 | Multi Clamp | Each | JP | φ16mm |

Side Rail Adaptor


Side Rail Adaptor, Long for Carbon Fiber Board



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|-----------------------|
| 08-113-01 | Side Rail Adaptor, Long for Carbon Fiber Board | Each | JP | |
| 08-113-01-R1 | Side Rail Adaptor, Long for Carbon Fiber Board (R1) | Each | US | |
| 08-113-01-R2 | Side Rail Adaptor, Long for Carbon Fiber Board (R2) | Each | EU | |


Side Rail Adaptor, Short for Carbon Fiber Board



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|---|
| 08-113-02 | Side Rail Adaptor, Short for Carbon Fiber Board | Each | JP | When attaching armboards, please use "Locking Type Armboard" on p.1 |
| 08-113-02-R1 | Side Rail Adaptor, Short for Carbon Fiber Board (R1) | Each | US | |
| 08-113-02-R2 | Side Rail Adaptor, Short for Carbon Fiber Board (R2) | Each | EU | |


Side Rail Adaptor II, Long for Carbon Fiber Board



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--|------|-----------|---|
| 08-113-10 | Side Rail Adaptor II, Long for Carbon Fiber Board | Each | JP | Knobs can be removed after attaching to the board |
| 08-113-10-R1 | Side Rail Adaptor II, Long for Carbon Fiber Board (R1) | Each | US | |
| 08-113-10-R2 | Side Rail Adaptor II, Long for Carbon Fiber Board (R2) | Each | EU | |

Side Rail Adaptor II, Short for Carbon Fiber Board




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|--|
| 08-113-11 | Side Rail Adaptor II, Short for Carbon Fiber Board | Each | JP | <ul style="list-style-type: none"> Knobs can be removed after attaching to the board When attaching armboards, please use "Locking Type Armboard" on p.1 |
| 08-113-11-R1 | Side Rail Adaptor II, Short for Carbon Fiber Board (R1) | Each | US | |
| 08-113-11-R2 | Side Rail Adaptor II, Short for Carbon Fiber Board (R2) | Each | EU | |

Headside Rail



- Auxiliary rail when the head section is removed
- Attach to the back section



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|------------------------------------|
| 08-113-50 | Headside Rail | Each | JP | Only hand controller is attachable |
| 08-113-50-R1 | Headside Rail (R1) | Each | US | |
| 08-113-50-R2 | Headside Rail (R2) | Each | EU | |

Protector for Hand Controller




- Non-slip silicon case provides extra grip
- Screen Protector provides added protection against external shock and damage

| Product Code | Description | Unit | Technical Information |
|--------------|--|------|--|
| 18-044-00 | Silicone Case for LCD Hand Controller | Each | Applicable table Type: MOT-VS700UKI |
| 18-044-10 | Screen Protector for LCD Hand Controller | Each | MST-7300 Series |

Instrument Tray



- Two joints allow for horizontal swivel
- Adjustable height




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|----------------------|------|-----------|---|
| 08-097-00-UE | Instrument Tray (UE) | Each | JP,US,EU | Requires Clamp for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Adjustable Instrument Tray



- Keeps the angle of tray when table is in Trendelenburg position
- Height adjustable using clamp




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|--|
| 08-098-01-NC | Adjustable Instrument Tray (NC) | Each | JP,US,EU | Requires two Clamps for Round Bar (Refer to p.45) JP : 08-110-01 US : 08-117-09 EU : 08-117-11 Universal : 08-110-20 |

Mayo Instrument Tray with Stand



- Flat legs are easy to insert under table base
- Adjustable height




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------------|------|-----------|-----------------------|
| 08-333-00 | Mayo Instrument Tray with Stand | Each | N/A | |

Accessory Cart (A)



- Multipurpose cart
- Easy transportation and storage of accessories



| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|--------------------|------|-----------|-----------------------|
| 08-081-01 | Accessory Cart (A) | Each | JP | |


Accessory Cart (B)



- Multipurpose cart
- Easy transportation and storage of table sections and accessories



*Accessories in the image not included




| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------------------|------|-----------|-----------------------|
| 08-081-09 | Accessory Cart (B) | Each | JP | |
| 08-081-09-R1 | Accessory Cart (B) (R1) | Each | US | |
| 08-081-09-R2 | Accessory Cart (B) (R2) | Each | EU | |

Table Cart



- Easy transportation and storage of table sections and larger accessories



(mm)




*Accessories in the image not included

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|-------------|------|-----------|-----------------------|
| 08-007-20 | Table Cart | Each | N/A | |


Cervical Spine Frame



- For prone procedures
- Radiolucent
- Use with spinal frame (refer to p.30 and p.31)




(mm)



(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---------------------------|------|-----------|--|
| 08-007-22 | Cervical Spine Frame | Each | JP | • Pad not included • Applicable tables: MOT-5701 MOT-VS700UKI MOT-VS600 Series MOT-VS500SK Series MST-7300 Series MOT-3602 *MOT-VS500DK Series *MOT-VE100D *When used with Leg-Side Reversal Adapter • Contents: Refer to ①-⑦ below |
| 08-007-22-R1 | Cervical Spine Frame (R1) | Each | US | |
| 08-007-22-R2 | Cervical Spine Frame (R2) | Each | EU | |


*Table attachment and head holder are not included. Only to be attached to Mizuho Head Holder.



Sitting Position Frame



- For use in neurosurgery procedures requiring sitting position
- Attached to table side rails




(mm)

| Product Code | Description | Unit | Side Rail | Technical Information |
|--------------|---|------|-----------|--|
| 07-953-01 | Sitting Position Frame with Side Rail Clamp | Each | JP | Requires two Clamps for Round Bar (Refer to p.45) US : 08-117-09 EU : 08-117-11 |

*Head holder is not included. Only to be attached to Mizuho Head Holder and Table.

SUSANO Table Pads



- 70-mm mattress made to flame-retardant* and waterproof** specifications
- Welded to keep out disinfectant and other liquids for improved hygiene and product durability
- Triple-layer mattress distributes body pressure more effectively
Upper/lower layers: soft urethane foam adapts to body contours
Middle layer: high-density foam supports body weight
- Latex-free

Choose required pads from item list (pp.54-58)

Note: Choose all sections for whole table configuration

* BS7177 Medium Hazard 5 certified

** IPX4 equivalent

Table Pads - Standard




- 70-mm latex-free mattress
- Triple-layer mattress distributes body pressure more effectively

Choose required pads from item list (pp.54-58)


Note: Choose all sections for whole table configuration

MOT-VS700UKI with Double Leg Section




| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-52 SUSANO Pad for Double Leg Section Type A | Pair | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-50 Mizuho Original Pad for Double Leg Section Type A | Pair | |

MOT-VS700UKI with Single Leg Section




| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-61 SUSANO Pad for Single Leg Section Type A | Each | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-60 Mizuho Original Pad for Single Leg Section Type A | Each | |

MOT-VS600 Series




| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-52 SUSANO Pad for Double Leg Section Type A | Pair | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-50 Mizuho Original Pad for Double Leg Section Type A | Pair | |

MOT-VS500DK




| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-53 SUSANO Pad for Double Leg Section Type B | Pair | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-51 Mizuho Original Pad for Double Leg Section Type B | Pair | |

MOT-VS500SK



| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-61 SUSANO Pad for Single Leg Section Type A | Each | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-60 Mizuho Original Pad for Single Leg Section Type A | Each | |

MOT-VE100D




| Product Code | Description | Unit | Technical Information |
|--|--|------|--|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-42 SUSANO Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-53 SUSANO Pad for Double Leg Section Type B | Pair | |
| Standard | | | |
| 1 | 18-042-00 Mizuho Original Pad for Head Section Type A | Each | |
| 2 | 18-042-40 Mizuho Original Pad for Back and Seat Section Type A | Each | Can be changed to split type I or II |
| 3 | 18-042-51 Mizuho Original Pad for Double Leg Section Type B | Pair | |

Back and Seat Section Split Type




| Product Code | Description | Unit | Technical Information |
|--|---|------|-----------------------|
| Susano - Flame retardant / Waterproof | | | |
| 2-1 | 18-042-11 SUSANO Pad for Back Section Type A | Each | |
| 2-2 | 18-042-22 SUSANO Pad for Seat Section Type A | Each | |
| Standard | | | |
| 2-1 | 18-042-10 Mizuho Original Pad for Back Section Type A | Each | |
| 2-2 | 18-042-20 Mizuho Original Pad for Seat Section Type A | Each | |

Back and Two Seat Section Split Type




| Product Code | Description | Unit | Technical Information |
|--|---|------|-----------------------|
| Susano - Flame retardant / Waterproof | | | |
| 2-1 | 18-042-11 SUSANO Pad for Back Section Type A | Each | |
| 2-2 | 18-042-23 SUSANO Pad for Seat Section Type B | Each | |
| Standard | | | |
| 2-1 | 18-042-10 Mizuho Original Pad for Back Section Type A | Each | |
| 2-2 | 18-042-21 Mizuho Original Pad for Seat Section Type B | Each | |

MOS-1302B




| Product Code | Description | Unit | Technical Information |
|--|--|------|-----------------------|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-13 SUSANO Pad for Back Section Type C | Each | |
| 3 | 18-042-25 SUSANO Pad for Seat Section Type D | Each | |
| 4 | 18-042-26 SUSANO Pad for Seat Section Type E | Pair | |
| 5 | 18-042-54 SUSANO Pad for Double Leg Section Type C | Pair | |
| Standard | | | |
| 08-055-60 | Pad for 1302B | Set | |

MST-7300 Series




| Product Code | Description | Unit | Technical Information |
|--|--|------|-----------------------|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-02 SUSANO Pad for Head Section Type B | Each | |
| 3 | 18-042-12 SUSANO Pad for Back Section Type B | Each | |
| 4 | 18-042-24 SUSANO Pad for Seat Section Type C | Each | |
| 5 | 18-042-62 SUSANO Pad for Single Leg Section Type B | Each | |
| Standard | | | |
| 08-055-70 | Pad for 7300 Series | Set | |

MOT-5602 Series




| Product Code | Description | Unit | Technical Information |
|--|---|------|-----------------------|
| Susano - Flame retardant / Waterproof | | | |
| 1 | 18-042-01 SUSANO Pad for Head Section Type A | Each | |
| 2 | 18-042-44 SUSANO Pad for Back and Seat Section Type C | Each | |
| 3 | 18-042-53 SUSANO Pad for Double Leg Section Type B | Pair | |
| Standard | | | |
| 08-055-25 | Pad for 5602BW & 5602SRMB | Set | |

VACB-3006B



| Product Code | Description | Unit | Technical Information |
|-----------------|---------------|------|-----------------------|
| Standard | | | |
| 08-055-75 | Pad for 3006B | Each | Thickness: 50mm |











SPL-331



| Product Code | Description | Unit | Technical Information |
|-----------------|-----------------|------|----------------------------------|
| Standard | | | |
| 08-010-10 | Pad for SPL-331 | Set | Thickness: 50mm *Single layer |

Power Cords

Cord can be selected according to each regional type.

| Product Code | Description | Unit | Example of Countries | Cord Length (m) | Image | Type | Standard |
|--------------|---------------|------|----------------------|-----------------|--|---|--|
| 18-045-02 | Power Cord 2 | Each | United States | 4.5 |  |  | ANSI/NEMA |
| 18-045-03 | Power Cord 3 | Each | Germany | 4.5 |  |  | VDE/SEMKO/ DEMKO/NEMKO/ FIMKO/CEBEC/ KEMA/OVE/CEE |
| 18-045-04 | Power Cord 4 | Each | Australia | 4.6 |  |  | AS/NZS |
| 18-045-05 | Power Cord 5 | Each | United Kingdom | 4.6 |  |  | BS |
| 18-045-06 | Power Cord 6 | Each | Denmark | 4.6 |  |  | DK |
| 18-045-07 | Power Cord 7 | Each | South Africa | 4.6 |  |  | SANS |
| 18-045-08 | Power Cord 8 | Each | Israel | 4.5 |  |  | SI |
| 18-045-09 | Power Cord 9 | Each | Argentina | 4.6 |  |  | IRAM |
| 18-045-10 | Power Cord 10 | Each | Italy | 4.6 |  |  | CEI |
| 18-045-11 | Power Cord 11 | Each | China | 4.5 |  |  | CCC |
| 18-045-12 | Power Cord 12 | Each | Switzerland | 4.6 |  |  | SEV |
| 18-045-13 | Power Cord 13 | Each | Philippines | 4.5 |  |  | NEMA |
| 18-045-14 | Power Cord 14 | Each | Brasil | 4.6 |  |  | NBR |

| Product Code | Description | Page No. |
|-----------------------|---|----------|
| 1. ARM SUPPORT | | |
| 08-065-70 | Locking Type Armboard | 1 |
| 08-065-80 | Locking Type Armboard (C1) | 1 |
| 08-065-85 | Locking Type Armboard (C2) | 1 |
| 08-065-08 | 40mm Pad for Armboard | 1 |
| 08-065-40 | 60mm Pad for Armboard | 1 |
| 08-065-45 | 100mm Pad for Armboard | 1 |
| 08-117-07 | 60mm Pad for Armboard (C1/C2) | 1 |
| 08-117-08 | 100mm Pad for Armboard (C1/C2) | 1 |
| 08-065-51 | Double Axis Armboard | 1 |
| 08-065-51-R1 | Double Axis Armboard (R1) | 1 |
| 08-065-51-R2 | Double Axis Armboard (R2) | 1 |
| 08-065-52 | 60mm Pad for Double Axis Armboard | 1 |
| 08-065-53 | 100mm Pad for Double Axis Armboard | 1 |
| 08-084-03 | Arm Strap | 2 |
| 08-065-72 | IVR Armboard | 2 |
| 08-066-06 | Raised Armboard with Clamp (730mm) | 3 |
| 08-066-10-UE | Raised Armboard (UE) | 3 |
| 08-066-05-NC | Armboard for Prone Position (NC) | 3 |
| 08-086-15 | Lower Armboard with Clamp | 4 |
| 08-086-15-R1 | Lower Armboard (R1) | 4 |
| 08-086-15-R2 | Lower Armboard (R2) | 4 |
| 08-065-13 | Imaging Armboard | 4 |
| 08-065-13-R1 | Imaging Armboard (R1) | 4 |
| 08-065-13-R2 | Imaging Armboard (R2) | 4 |
| 08-065-26 | 50mm Pad for Imaging Armboard | 4 |
| 08-065-27 | 90mm Pad for Imaging Armboard | 4 |
| 08-086-05 | Arm Surgery Board | 5 |
| 08-086-05-C1 | Arm Surgery Board with Side Rail Mount (C1) | 5 |
| 08-086-05-C2 | Arm Surgery Board with Side Rail Mount (C2) | 5 |
| 03-001-74 | Radiolucent Acrylic Armboard | 5 |
| 08-086-27 | Arm Shield with Pad (130mm) | 6 |
| 08-086-28 | Arm Shield with Pad (160mm) | 6 |
| 08-086-29 | Arm Shield with Pad (190mm) | 6 |
| 08-069-45 | Pad for Arm Shield (150mm) | 6 |
| 08-079-00 | Wrist Strap (240mm) | 6 |
| 08-079-00-R1 | Wrist Strap (240mm) (R1) | 6 |
| 08-079-00-R2 | Wrist Strap (240mm) (R2) | 6 |
| 08-068-04 | Shoulder Supports | 7 |
| 08-064-22-NC | Upper Arm Support II (NC) | 7 |
| 08-064-23-NC | Upper Arm Support Joint Type II (NC) | 8 |

| | | |
|-----------------------|--|---|
| 2. LEG SUPPORT | | |
| 18-070-91 | Levigator Stirrup II Set, Regular, Right Leg | 9 |
| 18-070-92 | Levigator Stirrup II Set, Regular, Left Leg | 9 |
| 18-070-93 | Levigator Stirrup II Main Unit, Right Leg | 9 |
| 18-070-94 | Levigator Stirrup II Main Unit, Left Leg | 9 |

Accessories List

| Product Code | Description | Page No. |
|--------------|--|----------|
| 18-070-50 | Short Boot Extender | 9 |
| 18-070-21 | Boot for Levitator Stirrup II, Regular, Right Leg | 10 |
| 18-070-22 | Boot for Levitator Stirrup II, Regular, Left Leg | 10 |
| 18-070-41 | Boot Pad for Levitator Stirrup II, Regular, Right Leg | 10 |
| 18-070-42 | Boot Pad for Levitator Stirrup II, Regular, Left Leg | 10 |
| 18-070-23 | Boot for Levitator Stirrup II, Fin Type, Right Leg | 10 |
| 18-070-24 | Boot for Levitator Stirrup II, Fin Type, Left Leg | 10 |
| 18-070-43 | Boot Pad for Levitator Stirrup II, Fin Type, Right Leg | 10 |
| 18-070-44 | Boot Pad for Levitator Stirrup II, Fin Type, Left Leg | 10 |
| 18-070-27 | Boot for Levitator Stirrup II, Wide, Right Leg | 10 |
| 18-070-28 | Boot for Levitator Stirrup II, Wide, Left Leg | 10 |
| 18-070-47 | Boot Pad for Levitator Stirrup II, Wide, Right Leg | 10 |
| 18-070-48 | Boot Pad for Levitator Stirrup II, Wide, Left Leg | 10 |
| 18-070-25 | Boot for Levitator Stirrup II, Short, Right Leg | 10 |
| 18-070-26 | Boot for Levitator Stirrup II, Short, Left Leg | 10 |
| 18-070-45 | Boot Pad for Levitator Stirrup II, Short, Right Leg | 10 |
| 18-070-46 | Boot Pad for Levitator Stirrup II, Short, Left Leg | 10 |
| 18-070-61 | Levitator Stirrups Cart | 11 |
| 08-070-57 | Transfer Board II (Rail Fixation Type) | 12 |
| 08-070-00-UE | Knee Crutches (UE) | 12 |
| 08-070-13 | Urethane Pads for Knee Crutches | 12 |
| 08-070-16-NC | Pediatric Knee Crutches, Large (NC) | 13 |
| 08-070-17-NC | Pediatric Knee Crutches, Middle (NC) | 13 |
| 08-070-18-NC | Pediatric Knee Crutches, Small (NC) | 13 |
| 08-072-01 | Knee Rests | 14 |
| 08-072-02 | Radiographic Tops for Knee Rests | 14 |
| 08-072-03 | 70mm Pads for Knee Rests | 14 |
| 08-071-03 | Foot Rests (320mm) | 14 |
| 08-071-05 | Mizuho Original Pads for Foot Rests | 14 |

3. HEAD SUPPORT

| | | |
|--------------|--------------------------------------|----|
| 08-073-03 | U-Shaped Head Rest, L | 15 |
| 08-073-04 | U-Shaped Head Rest, S | 15 |
| 07-952-00 | Table Attachment | 15 |
| 08-075-00 | Narrow Head Section | 16 |
| 08-075-00-R1 | Narrow Head Section (R1) | 16 |
| 08-075-00-R2 | Narrow Head Section (R2) | 16 |
| 08-083-69 | 70mm Pad for Narrow Head Section | 16 |
| 08-083-70 | 110mm Pad for Narrow Head Section | 16 |
| 08-075-06 | Head Supporter | 16 |
| 08-075-06-R1 | Head Supporter (R1) | 16 |
| 08-075-06-R2 | Head Supporter (R2) | 16 |
| 08-075-15 | Adapter for Carbon Fiber Extension | 17 |
| 08-075-11 | Attachment for Mayfield Head Support | 17 |

| Product Code | Description | Page No. |
|------------------------|------------------------------|----------|
| 4. BODY SUPPORT | | |
| 08-069-15-UE | Chest & Waist Support (UE) | 18 |
| 08-069-24-NC | Lateral Support (NC) | 18 |
| 08-080-11-NC | Adjustable Body Support (NC) | 19 |
| 08-084-02 | Body Strap (1700mm) | 19 |
| 08-084-00 | Leg Strap | 19 |
| 08-084-01 | Body Strap (2500mm) | 19 |

| | | |
|-----------------------------|---|----|
| 5. X-RAY ACCESSORIES | | |
| 18-090-00 | Radiographic Table Tops Fixed Type for Head Section Type A | 20 |
| 18-092-10 | Radiographic Table Tops Fixed Type for Back Section Type KA | 20 |
| 18-092-20 | Radiographic Table Tops Fixed Type for Seat Section Type KA | 20 |
| 18-090-50 | Radiographic Table Tops Fixed Type for Double Leg Section Type A | 20 |
| 18-090-60 | Radiographic Table Tops Fixed Type for Single Leg Section Type A | 20 |
| 18-090-10 | Radiographic Table Tops Fixed Type for Back Section Type A | 20 |
| 18-090-20 | Radiographic Table Tops Fixed Type for Seat Section Type A | 20 |
| 18-092-21 | Radiographic Table Tops Fixed Type for Seat Section Type KB | 21 |
| 18-090-51 | Radiographic Table Tops Fixed Type for Double Leg Section Type B | 21 |
| 18-090-61 | Radiographic Table Tops Fixed Type for Single Leg Section Type A-II | 21 |
| 08-089-72 | Radiographic Table Tops Fixed Type for 5602BW | 21 |
| 08-089-76 | Radiographic Table Tops Fixed Type for 7300 Series | 21 |
| 18-091-00 | Digital Radiographic Table Tops Fixed Type for Head Section Type A | 22 |
| 18-093-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type KA | 22 |
| 18-093-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KA | 22 |
| 18-091-60 | Digital Radiographic Table Tops Fixed Type for Single Leg Section Type A | 22 |
| 18-091-10 | Digital Radiographic Table Tops Fixed Type for Back Section Type A | 22 |
| 18-091-20 | Digital Radiographic Table Tops Fixed Type for Seat Section Type A | 22 |
| 18-093-21 | Digital Radiographic Table Tops Fixed Type for Seat Section Type KB | 23 |
| 18-091-61 | Digital Radiographic Table Tops Fixed Type for Single Leg Section Type A-II | 23 |
| 08-089-82 | Digital Radiographic Table Tops Fixed Type for 5602BW | 23 |
| 08-089-86 | Digital Radiographic Table Tops Fixed Type for 7300 Series | 23 |
| 08-088-20 | Radiographic Table Top for MOT-5602SRMB | 24 |
| 08-088-59 | Radiographic Table Tops for 1302B | 24 |
| 08-088-22 | Radiographic Table Tops for SPL-331 | 24 |
| 08-094-00 | Cassette Tray (390 x 466 mm) | 25 |
| 08-095-50 | Adjustable Cassette Tray | 25 |

| | | |
|--------------------------------------|---------------------------------------|----|
| 6. ORTHOPEDICS / TRAUMATOLOGY | | |
| 08-051-00 | Leg Traction Device MOC-142 | 26 |
| 08-115-01 | Traction Device Arm | 27 |
| 08-115-02 | Extension Arm, Right/Left | 27 |
| 08-115-04 | Traction Device Right/Left | 27 |
| 08-115-05 | Traction Shoe Right/Left | 27 |
| 08-115-07 | Inner Thigh Column with Pad | 27 |
| 08-115-27 | Support Rods | 27 |
| 08-115-76 | Wide Pelvic Support II | 27 |
| 08-115-77 | Transfer Board II (Pin Fixation Type) | 27 |

Accessories List

| Product Code | Description | Page No. |
|--------------|--|----------|
| 08-115-10 | Lower Limb Support | 27 |
| 08-115-11 | Inner Thigh Support with Pad for Traction in Lateral Posture | 27 |
| 08-115-12 | Lower Thigh Support for Traction in Lateral Posture | 27 |
| 08-115-26 | Pelvic Support Board | 27 |
| 08-115-26-R1 | Pelvic Support Board (R1) | 27 |
| 08-115-26-R2 | Pelvic Support Board (R2) | 27 |
| 08-115-09 | L-Shaped Knee Rest with Pad | 27 |
| 08-081-07 | Accessory Cart for Traction Device | 27 |
| 08-081-07-R1 | Accessory Cart for Traction Device (R1) | 27 |
| 08-081-07-R2 | Accessory Cart for Traction Device (R2) | 27 |
| 08-115-08 | Pelvic Support for Pediatric | 27 |
| 08-115-17 | Traction Pulley | 27 |
| 08-115-18 | Direct Traction Device | 27 |
| 08-115-19 | Direct Traction Hook Attachment | 27 |
| 08-115-20 | Movable Pelvic Support | 27 |
| 08-115-36 | Inner Thigh Column | 27 |
| 08-115-37 | Pad for Inner Thigh Column | 27 |
| 08-115-38 | Lower Limb Support Clamp | 27 |
| 08-115-43 | Carbon Fiber L Shaped Knee Rest with Pad | 27 |
| 08-115-45 | Direct Traction Device for Borshards | 27 |
| 08-115-73 | Lower Limb Support without Fixing Bracket | 27 |
| 08-115-21 | Cart for Traction Arm | 29 |
| 08-115-39 | Cart for Traction Arm Type II | 29 |
| 08-115-40 | Cart for Traction Arm Type III | 29 |
| 08-115-41 | Cart for Traction Arm Type IV | 29 |
| 08-115-42 | Cart for Traction Arm Type V | 29 |
| 08-115-55 | Cart for Traction Arm Type VI | 29 |
| 08-115-80 | V Attachment II Side Rail Type | 29 |
| 08-115-80-R1 | V Attachment II Side Rail Type (R1) | 29 |
| 08-115-80-R2 | V Attachment II Side Rail Type (R2) | 29 |
| 08-115-81 | V Attachment II Side Rail Type for Radiographic Table Top | 29 |
| 08-115-81-R1 | V Attachment II Side Rail Type Radiographic Table Top (R1) | 29 |
| 08-115-81-R2 | V Attachment II Side Rail Type Radiographic Table Top (R2) | 29 |
| 08-115-23 | V Attachment Pin Type | 30 |
| 08-115-23-R1 | V Attachment Pin Type (R1) | 30 |
| 08-115-23-R2 | V Attachment Pin Type (R2) | 30 |
| 08-115-24 | V Attachment Pin Type for Radiographic Table Top | 30 |
| 08-115-24-R1 | V Attachment Pin Type for Radiographic Table Top (R1) | 30 |
| 08-115-24-R2 | V Attachment Pin Type for Radiographic Table Top (R2) | 30 |
| 02-120-00 | GM-98 Spinal Frame | 30 |
| 02-121-00 | Four-Point Spinal Frame | 31 |
| 02-121-10 | Radiolucent Carbon Fiber Four Point Spinal Frame | 31 |
| 08-068-22 | Shoulder Positioner without Pad | 32 |
| 08-068-23 | Pad for Shoulder Positioner | 32 |
| 08-068-24 | Shoulder Positioner Replacement Strap for Head Support | 32 |
| 08-068-25 | Shoulder Positioner Replacement Pad, Head Section | 32 |
| 08-068-26 | Shoulder Positioner Replacement Pad, Shoulder Section | 32 |

| Product Code | Description | Page No. |
|--------------|---|----------|
| 08-068-27 | Shoulder Positioner Replacement Pad, Back Section | 32 |
| 08-078-32-NC | Pubis Support (NC) | 32 |

7. EXTENSION / SUPPORT ACCESSORIES

| | | |
|--------------|---|----|
| 08-086-06 | Table Side Extension | 33 |
| 08-086-06-R1 | Table Side Extension (R1) | 33 |
| 08-086-06-R2 | Table Side Extension (R2) | 33 |
| 08-086-33 | 70mm Pad for Table Side Extension | 33 |
| 08-086-35 | 110mm Pad for Table Side Extension | 33 |
| 08-075-90 | Head Extension II | 33 |
| 08-075-90-R1 | Head Extension II (R1) | 33 |
| 08-075-90-R2 | Head Extension II (R2) | 33 |
| 08-083-02 | Radiographic Table Tops Fixed Type for Head Extension (240mm) | 33 |
| 08-083-22 | Digital Radiographic Table Tops Fixed Type for Head Extension (240mm) | 33 |
| 18-042-01 | SUSANO Pad for Head Section Type A | 33 |
| 18-042-00 | Mizuho Original Pad for Head Section Type A | 33 |
| 08-075-36 | Seat Extension (350mm) | 34 |
| 08-075-36-R1 | Seat Extension (350mm) (R1) | 34 |
| 08-075-36-R2 | Seat Extension (350mm) (R2) | 34 |
| 08-075-42 | Seat Extension for 5701 | 34 |
| 08-075-42-R1 | Seat Extension for 5701 (R1) | 34 |
| 08-075-42-R2 | Seat Extension for 5701 (R2) | 34 |
| 08-083-01 | Radiographic Table Tops Fixed Type for Seat Extension (350mm) | 34 |
| 08-083-21 | Digital Radiographic Table Tops Fixed Type for Seat Extension (350mm) | 34 |
| 08-083-00 | Radiographic Table Tops Fixed Type for Seat Extension (300mm) | 34 |
| 08-083-20 | Digital Radiographic Table Tops Fixed Type for Seat Extension (300mm) | 34 |
| 18-042-83 | SUSANO Pad for Seat Extension (350mm) | 34 |
| 18-042-82 | SUSANO Pad for Seat Extension (300mm) | 34 |
| 08-083-61 | Mizuho Original Pad for Seat Extension (350mm) | 34 |
| 08-083-60 | Mizuho Original Pad for Seat Extension (300mm) | 34 |
| 08-075-52 | Seat Extension II (300mm) | 35 |
| 08-075-52-R1 | Seat Extension II (300mm) (R1) | 35 |
| 08-075-52-R2 | Seat Extension II (300mm) (R2) | 35 |
| 08-075-30 | Leg Extension for 5700 Series | 36 |
| 08-075-30-R1 | Leg Extension for 5700 Series (R1) | 36 |
| 08-075-30-R2 | Leg Extension for 5700 Series (R2) | 36 |
| 08-083-04 | Radiographic Table Tops Fixed Type for Leg Extension | 36 |
| 18-042-84 | SUSANO Pad for Leg Extension | 36 |
| 08-083-62 | Mizuho Original Pad for Leg Extension | 36 |
| 08-115-32 | Side Rail Extension, 220mm (Pair) | 37 |
| 08-095-08 | Uro Drain Tray Set for 5602SRMB | 37 |
| 08-095-13 | Uro Drain Tray Set for 3600 & 3500 Series, 5801 | 37 |
| 08-095-15 | Uro Drain Tray Set for 5701 | 37 |
| 08-095-19 | Uro Drain Tray for SPL-331 | 37 |
| 08-095-20 | Waste Container with Hose | 37 |
| 08-068-31 | Back Section Board Attachment for Pediatric | 38 |
| 08-068-31-R1 | Back Section Board Attachment for Pediatric (R1) | 38 |

Accessories List

| Product Code | Description | Page No. |
|--------------|--|----------|
| 08-068-31-R2 | Back Section Board Attachment for Pediatric (R2) | 38 |
| 08-083-06 | Radiographic Table Tops Fixed Type for Back Section Board Attachment for Pediatric | 38 |
| 18-042-81 | SUSANO Pad for Back Section Board Attachment for Pediatric | 38 |
| 08-083-63 | Mizuho Original Pad for Back Section Board Attachment for Pediatric | 38 |
| 08-007-21 | Back Section Board Short | 38 |
| 08-007-21-R1 | Back Section Board Short (R1) | 38 |
| 08-007-21-R2 | Back Section Board Short (R2) | 38 |
| 08-083-03 | Radiographic Table Tops Fixed Type for Back Section Board Short | 38 |
| 08-083-23 | Digital Radiographic Table Tops Fixed Type for Back Section Board Short | 38 |
| 18-042-80 | SUSANO Pad for Back Section Board Short | 38 |
| 08-083-64 | Mizuho Original Pad (252mm) | 38 |
| 08-096-02 | Kidney Bridge with Pad | 39 |
| 08-070-54 | 100cm Carbon Fiber Extension II | 39 |
| 08-083-65 | 70mm Mizuho Original Pad for 100cm Carbon Fiber Extension | 39 |
| 08-083-66 | 110mm Mizuho Original Pad for 100cm Carbon Fiber Extension | 39 |
| 08-070-55 | 122cm Carbon Fiber Extension II | 40 |
| 08-083-67 | 70mm Mizuho Original Pad for 122cm Carbon Fiber Extension | 40 |
| 08-083-68 | 110mm Mizuho Original Pad for 122cm Carbon Fiber Extension | 40 |
| 08-077-70 | Leg-Side Reversal Adapter A | 41 |
| 08-077-60 | Leg-Side Reversal Adapter B | 41 |
| 08-077-60-R1 | Leg-Side Reversal Adapter B (R1) | 41 |
| 08-077-60-R2 | Leg-Side Reversal Adapter B (R2) | 41 |
| 08-077-50 | Leg-Side Reversal Adapter C | 41 |
| 08-077-50-R1 | Leg-Side Reversal Adapter C (R1) | 41 |
| 08-077-50-R2 | Leg-Side Reversal Adapter C (R2) | 41 |
| 18-097-00 | Digital Radiographic Table Top for Leg-Side Reversal Adapter | 41 |
| 18-042-90 | Mizuho Original Pad for Leg-Side Reversal Adapter | 41 |
| 18-443-43 | Foot Switch | 42 |
| 18-443-50 | Foot Switch Type A | 42 |
| 18-443-20 | Foot Switch Type B | 42 |

8. GENERAL ACCESSORIES

| | | |
|--------------|---|----|
| 08-062-00 | Anesthesia Screen L Type with Clamp | 43 |
| 08-062-00-R1 | Anesthesia Screen L Type with Clamp (R1) | 43 |
| 08-062-00-R2 | Anesthesia Screen L Type with Clamp (R2) | 43 |
| 08-061-40 | Anesthesia Screen Extension with Clamp | 43 |
| 08-061-40-R1 | Anesthesia Screen Extension with Clamp (R1) | 43 |
| 08-061-40-R2 | Anesthesia Screen Extension with Clamp (R2) | 43 |
| 08-061-20-UE | Anesthesia Screen Wing Type (UE) | 44 |
| 08-110-01 | Clamp for Round Bar | 45 |
| 08-117-09 | Clamp for Round Bar (R1) | 45 |
| 08-117-11 | Clamp for Round Bar (R2) | 45 |
| 09-228-12 | Side Rail Clamp (φ18mm) | 45 |
| 08-110-20 | Universal Clamp for Round Bar | 45 |
| 08-111-11 | Clamp for Flat Bar Ver.2 | 45 |
| 08-117-10 | Clamp for Flat Bar Ver.2 (R1) | 45 |
| 08-117-13 | Clamp for Flat Bar Ver.2 (R2) | 45 |

| Product Code | Description | Page No. |
|--------------|---|----------|
| 08-111-20 | Universal Clamp for Flat Bar | 45 |
| 08-110-11 | Multi Clamp | 45 |
| 08-113-01 | Side Rail Adaptor, Long for Carbon Fiber Board | 46 |
| 08-113-01-R1 | Side Rail Adaptor, Long for Carbon Fiber Board (R1) | 46 |
| 08-113-01-R2 | Side Rail Adaptor, Long for Carbon Fiber Board (R2) | 46 |
| 08-113-02 | Side Rail Adaptor, Short for Carbon Fiber Board | 46 |
| 08-113-02-R1 | Side Rail Adaptor, Short for Carbon Fiber Board (R1) | 46 |
| 08-113-02-R2 | Side Rail Adaptor, Short for Carbon Fiber Board (R2) | 46 |
| 08-113-10 | Side Rail Adaptor II, Long for Carbon Fiber Board | 46 |
| 08-113-10-R1 | Side Rail Adaptor II, Long for Carbon Fiber Board (R1) | 46 |
| 08-113-10-R2 | Side Rail Adaptor II, Long for Carbon Fiber Board (R2) | 46 |
| 08-113-11 | Side Rail Adaptor II, Short for Carbon Fiber Board | 46 |
| 08-113-11-R1 | Side Rail Adaptor II, Short for Carbon Fiber Board (R1) | 46 |
| 08-113-11-R2 | Side Rail Adaptor II, Short for Carbon Fiber Board (R2) | 46 |
| 08-113-50 | Headside Rail | 47 |
| 08-113-50-R1 | Headside Rail (R1) | 47 |
| 08-113-50-R2 | Headside Rail (R2) | 47 |
| 18-044-00 | Silicone Case for LCD Hand Controller | 47 |
| 18-044-10 | Screen Protector for LCD Hand Controller | 47 |

9. TRAYS

| | | |
|--------------|---------------------------------|----|
| 08-097-00-UE | Instrument Tray (UE) | 48 |
| 08-098-01-NC | Adjustable Instrument Tray (NC) | 48 |
| 08-333-00 | Mayo Instrument Tray with Stand | 49 |

10. CARTS

| | | |
|--------------|-------------------------|----|
| 08-081-01 | Accessory Cart (A) | 50 |
| 08-081-09 | Accessory Cart (B) | 50 |
| 08-081-09-R1 | Accessory Cart (B) (R1) | 50 |
| 08-081-09-R2 | Accessory Cart (B) (R2) | 50 |
| 08-007-20 | Table Cart | 51 |

11. NEURO

| | | |
|--------------|---|----|
| 08-007-22 | Cervical Spine Frame | 52 |
| 08-007-22-R1 | Cervical Spine Frame (R1) | 52 |
| 08-007-22-R2 | Cervical Spine Frame (R2) | 52 |
| 07-953-01 | Sitting Position Frame with Side Rail Clamp | 52 |

12. TABLE PADS

| | | |
|-----------|--|----|
| 18-042-42 | SUSANO Pad for Back and Seat Section Type A | 54 |
| 18-042-52 | SUSANO Pad for Double Leg Section Type A | 54 |
| 18-042-40 | Mizuho Original Pad for Back and Seat Section Type A | 54 |
| 18-042-50 | Mizuho Original Pad for Double Leg Section Type A | 54 |
| 18-042-61 | SUSANO Pad for Single Leg Section Type A | 54 |
| 18-042-60 | Mizuho Original Pad for Single Leg Section Type A | 54 |
| 18-042-53 | SUSANO Pad for Double Leg Section Type B | 55 |
| 18-042-51 | Mizuho Original Pad for Double Leg Section Type B | 55 |

Accessories List

| Product Code | Description | Page No. |
|--------------|---|----------|
| 18-042-11 | SUSANO Pad for Back Section Type A | 56 |
| 18-042-22 | SUSANO Pad for Seat Section Type A | 56 |
| 18-042-10 | Mizuho Original Pad for Back Section Type A | 56 |
| 18-042-20 | Mizuho Original Pad for Seat Section Type A | 56 |
| 18-042-23 | SUSANO Pad for Seat Section Type B | 56 |
| 18-042-21 | Mizuho Original Pad for Seat Section Type B | 56 |
| 18-042-13 | SUSANO Pad for Back Section Type C | 57 |
| 18-042-25 | SUSANO Pad for Seat Section Type D | 57 |
| 18-042-26 | SUSANO Pad for Seat Section Type E | 57 |
| 18-042-54 | SUSANO Pad for Double Leg Section Type C | 57 |
| 08-055-60 | Pad for 1302B | 57 |
| 18-042-02 | SUSANO Pad for Head Section Type B | 57 |
| 18-042-12 | SUSANO Pad for Back Section Type B | 57 |
| 18-042-24 | SUSANO Pad for Seat Section Type C | 57 |
| 18-042-62 | SUSANO Pad for Single Leg Section Type B | 57 |
| 08-055-70 | Pad for 7300 Series | 57 |
| 18-042-44 | SUSANO Pad for Back and Seat Section Type C | 58 |
| 08-055-25 | Pad for 5602BW & 5602SRMB | 58 |
| 08-055-75 | Pad for 3006B | 58 |
| 08-010-10 | Pad for SPL-331 | 58 |

13. POWER CORDS

| | | |
|-----------|---------------|----|
| 18-045-02 | Power Cord 2 | 59 |
| 18-045-03 | Power Cord 3 | 59 |
| 18-045-04 | Power Cord 4 | 59 |
| 18-045-05 | Power Cord 5 | 59 |
| 18-045-06 | Power Cord 6 | 59 |
| 18-045-07 | Power Cord 7 | 59 |
| 18-045-08 | Power Cord 8 | 59 |
| 18-045-09 | Power Cord 9 | 59 |
| 18-045-10 | Power Cord 10 | 59 |
| 18-045-11 | Power Cord 11 | 59 |
| 18-045-12 | Power Cord 12 | 59 |
| 18-045-13 | Power Cord 13 | 59 |
| 18-045-14 | Power Cord 14 | 59 |



MIZUHO Corporation

3-30-13 Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN

Phone +81-3-6627-1506

F a x +81-3-3815-1280

U R L <https://www.mizuho.co.jp/en>



Website



Contact Us

* Under our policy, ongoing improvements and published specifications may be subject to change without prior notice.

* All published CG photos and illustrations have been from images.

©MIZUHO 2025



EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands



| | |
|--|---|
|  Caution | <ul style="list-style-type: none">• Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.• Inspect the product routinely and thoroughly, in addition to pre- and postoperative examinations.• For after-sale support, contact your local dealer / distributor. |
|--|---|

